

CleanAmp™ Hot-Start 7-deaza-dGTP for Improved GC-rich PCR Amplification

Sabrina Shore, Natasha Paul, Richard Hogrefe, Gerald Zon, Alexandre Lebedev

Abstract

PCR amplification of nucleic acids is a fundamental technique used in many molecular biology laboratories. Despite its widespread use, GC-rich regions of DNA sequence remain a challenge for amplification. Sequences high in GC content can form strong secondary structure, which prevents strand denaturation and blocks processive DNA polymerase amplification. As a consequence, mis-priming is prominent and can complicate specific product formation. Especially as applied to the molecular diagnosis of inheritable diseases, several assay modifications have been developed to improve the specificity of target amplification. These approaches include specialized polymerases, Hot Start assays, addition of organic molecules and thermal cycling alterations. However, as the GC content increases, the combination of two or three approaches may be required. Here, we show how 7-deaza-dGTP, a commonly used molecule to amplify GC-rich targets, can greatly improve results when a thermolabile protecting group is incorporated at the 3'-hydroxyl. The presence of the protecting group blocks low temperature primer extension and only allows nucleotide incorporation at higher temperatures when the protecting group is removed, improving PCR specificity as a result. This Hot Start version of 7-deaza-dGTP, CleanAmp™ 7-deaza-dGTP improves the amplification of targets containing up to 80% GC content. Results were further improved when a Hot Start version of all dNTPs was employed, allowing for challenging targets of more than 85% GC content, such as Fragile X, to be amplified. Another benefit of this technology is in downstream sequencing reactions. PCR amplification of problematic targets with a Hot Start 7-deaza-dGTP mix prior to Sanger dideoxy sequencing can significantly improve the read quality along the entire sequence. In summary, the use of CleanAmp™ dNTPs simplifies GC-rich amplification and provides a valuable solution that can improve disease diagnosis.

Figure 1: CleanAmp™ Hot Start dNTP Chemical Structure




Figure 2: CleanAmp™ dNTP Activation Mechanism




Figure 3: 7-deaza-dGTP Reduces Secondary Structure Formation




Figure 4: Improved Amplification Yield and Specificity for Targets Up to 85% GC-Content




Figure 5: CleanAmp™ 7-deaza-dGTP Outperforms Common Additives for GC-Rich Target Amplification




Figure 6: Increased Specificity and Amplicon Yield Across a Range of Template Concentrations




Figure 7: Successful Multiplex PCR of GC-Rich Targets




Figure 8: Bacterial Identification by Amplification of Variable Number of Tandem Repeats (VNTRs)




Figure 9: PCR Amplification of GC-Rich Targets Improves the Accuracy of Dideoxy Sequencing




Figure 10: PCR Amplification Also Improves Sequencing Accuracy of Homopolymer Regions





Figure 11: PCR Amplification Allows for Longer Dideoxy Sequencing Results



Conclusion

- CleanAmp™ 7-deaza-dGTP Mix allows for clean amplification of high GC-rich targets of up to 85% GC content by reducing off-target amplification and increasing amplicon yield.
- CleanAmp™ 7-deaza-dGTP Mix outperforms other GC-rich amplification solutions.
- CleanAmp™ 7-deaza-dGTP Mix improves amplification efficiency and limit of detection of targets with high GC content while maintaining a clean and robust product.
- CleanAmp™ 7-deaza-dGTP Mix is able to successfully amplify three GC-rich targets in multiplex PCR.
- A PCR step with CleanAmp™ 7-deaza-dGTP Mix prior to dideoxy sequencing improves the quality of sequencing data for targets with high GC content by reducing background and improving base-calling.

Acknowledgements

Shawna Rumrill, Dave Combs, Terry Beck, Olga Adelfinskaya, Jessica Stewart, Elizabeth Hill, Angela Tenenini, Anton McCaffrey, Julie Powers.

NIH grants:

1R43GM079836-01, 2R44GM079836-02 and 5R44GM079836-03

Contact: Natasha Paul, npaul@trilinkbiotech.com
trilinkbiotech.com/cleanamp

LEGEND:

- 1X Standard 7-deaza-dGTP mix: 0.2 mM d(A, C, T, T)TP; 0.05 mM dGTP and 0.15 mM 7-deaza-dGTP
- 1X CleanAmp™ 7-deaza-dGTP Mix: 0.2 mM CleanAmp™ d(A, C, T, T)TP; 0.05 mM CleanAmp™ dGTP and 0.15 mM CleanAmp™ 7-deaza-dGTP