## Novel N1-Substituted Pseudouridine 5'-Triphosphates for the Synthesis of Modified mRNA and its Effect on mRNA Translation in THP-1 Cells



Dongwon Shin, Krist T. Azizian, Sabrina Shore, Jordana M. Henderson, Anton P. McCaffrey, Richard I. Hogrefe and Alexandre Lebedev TriLink BioTechnologies, San Diego, CA 92121, USA.

## Abstract

Messenger RNA (mRNA) therapeutics is an emerging platform technology for the treatment of gene disorders, immunotherapy and vaccines. It was that fested in 1990, but ubusequently found nonutable because of thin innel immunogenicity and to welfakercy. None executivy, Karlos et al. demonstrated that mRNA exonationing naturally occurring pseudourdine (%) increased overall protein agreesion advisorial control or protein strategies (RVS) and RNE-messide overall protein segression advisorial or protein strategies (RVS) and RNE-messide phosphore) advisorial or transition initiation factor 2 dapta (eR2a), Alchiedan of innate immune pathways leads to had dawn o potein transition Andres *et al.* channel that NH mRNysequecuration (RNA upon transfection into cell lines o mice.

electronic properties of the NI-substitution group of NI-substituted VPIA. The Fluck \*mRNMs with NI-substitution were fetted for expression of functionaries in the THF-1 monocyte cell line, a samitilite model for innote immune actinution. Four of seven NI-substituted mRNAs showed columbies higher than the activity of NI-substitutied VMRNA and ware close to the activity of NIsouther states that the activity of NI-substitutied VMRNA and ware close to the activity of NIsouther states that the activity of NI-substitutied VMRNA and ware close the activity of NIsouther states that the activity of NI-substitutied VMRNA and ware close the activity of NIsouther gene metacrostic that the state of the NI-substitutied VMRNA and ware close that where gene metacrostic throation states were the states that the states were active to the NI-substitutied where gene metacrostic throation activity ware due not activity the NI-substitutied activity and the states activity to THRI-substitutied activity activity activity activity activity activity the state of activity the NI-substitutied activity activity activity activity activity activity the state of activity the state of activity the NI-substitutied activity the state of activity the NI-substitutied activity the NI-substituties activity acti

#### • Background: Why mRNA therapeutics?

- mRNA is a popular new tool for gene expression because it:
- Does not have a risk of insertional mutagenesis
  Can transfect difficult cells such as non-dividing cells
- Is transient

#### Applications

- Genome editing (Transposons, Cre, ZFNs, TALENs and CRISPR/ Cas9)
- Gene replacement - Vaccines
- Limitations
- Innate immune response to unmodified mRNA

#### Solutions

- Proper capping
- Chemical modification of mRNA can prevent innate immune stimulation - Removal of dsRNA

# Innate immune sensors recognize mRNA

- Transfection of cells with unmodified RNAs can lead to cell death due to activation of innate immune pathways
- Toll-like receptors 3, 7 & 8 recognize different RNA forms
- Found in endosomes where some viruses enter cells

#### Cytosolic sensors

- Protein Kinase R (PKR): dsRNA - MDA5: long dsRNA - IFITs: unmethylated cap structures - RIG-I: 5'-triphosphate



- Eukaryotic mRNAs have a Cap1 or Cap2 structure.
- Sensing of proper cap structure is thought to be involved in self/ non-self RNA recognition.



 Co-transcriptional capping with CleanCap™ (Cap1) helps evade an immune response



Figure 3: Pseudouridine 5'-Triphosphate Derivatives

- Incorporation of modified nucleosides in mRNA helps to evade
  an immune response
- $\bullet$  Several novel pseudouridine NTPs were synthesized and tested in firefly luciferase mRNA transcriptions  $$\rm O$$



Pseudouridine 5-triphosphate derivatives; H = pseudouridine, Me = N1-methyl, EI = N1-ethyl, FE = fluoroethyl, P = sicopropt, MCM = methoxy methyl, POM = pivaloxy methyl, BOM = berzyloxy methyl pseudouridine



Figure 5: U Depletion of Primary Luciferase mRNA Sequence Improves Incorporation of N1-Substituted Ψ Derivatives by T7 Polymerase

- + Some N1-substituted  $\Psi$  derivatives did not incorporate well in WT mRNA
- We decreased the number of uridine residues in the sequence by substituting synonymous codons
- $\boldsymbol{\cdot}$  U depletion resulted in good incorporation



Figure 6: Sequence Engineering of FLuc mRNA (Ψ

or MOM1<sup>4</sup> Substitution, Bioanalyzer)



Figure 7: In Vitro Translation and Cell Activity of Modified Luciferase mRNAs



Figure 8: N1-Substituted  $\Psi$  Derivatives Resulted in Lower Toxicity Compared to WT and  $\Psi$  in FLuc mRNAs



### Conclusions

- We have synthesized a number of 5'-triphosphates of N1-modified pseudouridine derivatives
- These nucleoside 5'-triphosphates were used for the synthesis of modified mRNAs by *in vitro* transcription using WT and U-depleted templates
- Efficiency of transcription using U-depleted templates greatly improved quality and yield
- + N1-substituted  $\Psi$  mRNAs show potential translational and immunological properties
- Translational activity of modified mRNAs in wheat germ extracts did not directly correlate with cell activity, which may indicate differences in immune stimulation by these mRNAs

### Contact

Dongwon Shin dshin@trilinkbiotech.com

The Modified Nucleic SCI Acid Experts® Cze www.trilinkbiotech.com

SCNAC 2017 Czech Republic 06/17