

Maximizing Translation of mRNA Therapeutics By Sequence Engineering and Chemical Modification

Anton P. McCaffrey¹, Krist T. Azizian¹, Dongwon Shin¹, Sabrina Shore¹, Jordana M. Henderson¹, Alexandre Lebedev¹, Joel Jesse², Richard I. Hogrefe¹, Gerald Zon¹

¹TriLink BioTechnologies, San Diego, CA 92121, USA. ²MTI-GlobalStem, Gaithersburg MD 20877, USA.

Abstract

For maximal expression in cells or target organs, transfected mRNAs must avoid detection by pattern recognition receptors (PRRs) that evolved to sense pathogenic non-self RNAs. These include PRRs that recognize improperly capped RNAs (RG4, IRT1) and double stranded RNA (PKR, OAS, RIG-I, MDA5). PRR activation leads to cytokine production, translational arrest and cell toxicity or death. Mammalian mRNAs are modified post-transcriptionally to contain nucleotides with 2'-O-methyl residues, pseudouridine (Ψ) and N⁶-methyladenosine (m⁶A). Interestingly, these modifications can reduce activation of PRRs and allow maximal translation of the transfected mRNA. For mRNA drugs to achieve their potential, methods to produce multigram batches of properly capped, highly purified and non-immunogenic mRNA will be needed. For vaccines, some immunogenicity may be desirable to serve as an adjuvant.

During RNA capping, Cap0 (m⁷GpppN) is formed as an intermediate. Methylation of the 2' position of the first and sometimes second nucleotide forms Cap1 and Cap2, respectively. Recombinant enzymes used to generate Cap1 mRNA are expensive, do not always go to completion and the RNA must be purified prior to capping. We developed a novel co-transcriptional capping method called CleanCap™ that yields Cap1 with high efficiency and lower costs in a "one pot" reaction. CleanCap also facilitates modification of the 5' end of mRNAs with diverse functional groups, including m⁶A. We developed a capping assay that allows direct assessment of mRNA capping. Capping efficiencies as high as 99% can be obtained.

Previously, we identified 5-methyluridine (5mU) as a promising modification to avoid innate immune stimulation while supporting efficient translation. To further explore chemically modified bases, we synthesized several Ψ derivatives: N1-ethyl-Ψ (E1Ψ), N1-fluoroethyl-Ψ (FE1Ψ), N1-propyl-Ψ (P1Ψ), N1-isopropyl-Ψ (IP1Ψ) and N1-methyl-Ψ (MOM1Ψ). 5'-triphosphates (Luciferase mRNAs were fully substituted with these modifications and translational potential was monitored in wheat germ extracts. Activity was also measured in the THP-1 monocytic cell line, which is a sensitive model for innate immune activation. A recent report showed that minimizing uridine content in mRNAs reduced innate immune stimulation by unmethylated mRNAs. Here we show that incorporation of our modified uridine residues or 5mU into uridine depleted luciferase resulted in equivalent activity relative to wt THP-1 cells. We are currently evaluating the use of 5mU in depleted Renilla, beta-galactosidase, erythropoietin, mCherry and Cas9 mRNAs.

Background: Why mRNA therapeutics?

- mRNA is a popular new tool for gene expression because it:
 - Does not have a risk of insertional mutagenesis
 - Can transfect difficult cell types such as non-dividing cells
 - Is transient

Applications

- Genome editing (Transposons, Cre, ZFNs, TALENs and CRISPR/Cas9)
- Gene replacement
- Vaccines

Limitations

- Innate immune response to unmethylated mRNA

Solutions

- Proper capping
- Chemical modification of mRNA can prevent innate immune stimulation
- Removal of dsRNA

Innate immune sensors recognize mRNA

- Transfection of cells with unmethylated RNAs can lead to cell death due to activation of innate immune pathways

- Toll-like receptors 3, 7 & 8 recognize different RNA forms
 - Found in endosomes where some viruses enter cells

Cytosolic sensors

- Protein Kinase R (PKR): dsRNA
- MDA5: dsRNA
- IRT1: unmethylated cap structures
- RIG-I: 5' triphosphate

Figure 1: Cap0, Cap1 and Cap2 Structures of 5'-Ends of mRNAs

- Eukaryotic mRNAs have a Cap1 or Cap2 structure.
- Sensing of proper cap structure is thought to be involved in self/non-self RNA recognition.
- Cap structure influences activation of PRRs
 - RIG-I is activated by Cap0 RNAs but not Cap1 mRNAs (PMID: 18426922 and 20457754)
 - IRT1 binds Cap0 RNAs more tightly than Cap1 mRNAs (PMID: 24371270)






Figure 2: Capping Efficiency Assay Shows CleanCap™ Yields High Levels of Cap1



Capping efficiencies are consistently above 90%, depending on reaction conditions. Data collected by liquid chromatography with mass confirmed by mass spectrometry assay

Figure 3: Pseudouridine 5'-Triphosphate Derivatives

- mRNA body modifications help to evade an immune response
- Pseudouridine or 5-methylcytidine/pseudouridine are current industry standard
- Several novel pseudouridine NTPs were synthesized and tested in firefly luciferase transcriptions.



Pseudouridine 5'-triphosphate derivatives: H = pseudouridine, Me = N1-Methyl, B = N1-Ethyl, FE = Fluoroethyl, P = Propyl, IP = Isopropyl, MOM = Methylomethyl, POM = Pivalomethyl, BOM = Benzyloxymethyl pseudouridine

Figure 4: U Depletion of Primary Luciferase Sequence Improves Incorporation of Bulky Pseudouridine Derivatives by T7 Polymerase

- Some pseudouridine derivatives did not incorporate well
- We depleted the Fluc sequence for Us to try and remedy this
- U depletion resulted in good incorporation
- We tested the derivatives that did incorporate for translation and activity





Figure 5: In Vitro Translation and Cell Activity of Modified Luciferase mRNAs

- U depleted sequences translated better in wheat germ extracts
- Bulky pseudouridine modifications did not translate well
- U depleted sequences resulted in higher activity in THP-1 cells
- We therefore continued our studies using the U depleted sequence




Figure 6: Pseudouridine Derivatives and 5mU Resulted in Lower Toxicity Compared to WT and PseudoU




Figure 7: Cell Activity of GFP Modified RNAs in Various Cell Lines




Figure 8: Cell Activity of Standard and U Depleted Firefly Luciferase in THP-1 Dual Cells




Figure 9: Effect of HPLC Purification on Firefly Luciferase Cell Activity in THP-1 Dual Cells




Figure 10: Cell Activity of U Depleted Renilla Luciferase in THP-1 Dual Cells




Figure 11: Interferon Reporter Activity of THP-1 Dual Cells in Response to RNA Transfection





Figure 12: Slot Blot Demonstrates that HPLC Purification Depletes dsRNA

- An RP-HPLC method depletes mRNAs of contaminating dsRNA
- This reduces the innate immune response by reducing PKR activation



- As previously described, HPLC purification depleted dsRNA, reduced toxicity and reduced interferon activation





Figure 13: Cell Activity of HPLC vs. non-HPLC Luciferase mRNAs

- HPLC purification dramatically increased the activity of wt mRNA, improved the activity of PsU mRNA but did not alter the activity of 5mU mRNA
- Could this be because PKR does not bind 5mU mRNA?
- Could this also be true for the PseudoU derivatives?



Conclusions

- We have introduced a number of novel modified bases with interesting translational and immunological properties
- U depletion improved transcription quality, yield and activity
- HPLC purification to remove dsRNA reduced toxicity and interferon response and increased activity
- Translational activity in wheat germ extracts did not directly correlate with cell activity, which may indicate differences in immune stimulation by these mRNAs
- 5mU is a promising modification for reducing innate immune stimulation
- Ability of 5mU to suppress innate immune stimulation is sequence context dependent
- HPLC improves activity of WT and PseudoU modified RNAs but not 5mU modified RNAs. HPLC may not be necessary for 5mU
- One possibility is that 5mU dsRNA is not efficiently recognized by PRRs

Future Directions

- If 5mU is not recognized by PRRs, then activity of 5mU should be equivalent in PKR +/- and +/- MEFs
- Measure activity, toxicity and interferon response in THP-1 cells for HPLC purified PseudoU derivative mRNAs

Contact

Gerald Zon gzon@trilinkbiotech.com