CleanCap[®] Co-transcriptional Capping Streamlines mRNA Manufacturing

Dongwon Shin¹, Krist T. Azizian², Jordana M. Henderson¹, Richard I. Hogrefe¹, Alexandre Lebedev¹, Michael Houston¹, and Anton P. McCaffrey¹

Abstract

Messenger RNA (mRNA) therapy is a popular platform technology for expressing proteins in cells or *in vivo* because there is minimal risk of insertional mutagenesis. mRNA transfection is used to express proteins for genome editing, protein replacement, vaccines and antibody expression. To avoid an innate immune response, transfected mRNAs should mimic the 5' cap structure of non-immunogenic endogenous mRNAs.

During eukaryotic RNA capping, Cap 0 (^{m7}GpppN) is formed as an intermediate. Methylation of the 2'-O position of the first cap-proximal nucleotide forms Cap 1 (m7 GpppN $_{m}$ N). In ~50% of transcripts, the 2'-O position of the second cap-proximal nucleotide is also methylated to form Cap 2 (^{m7}GpppN_mN_m). N6-methylation of adenosine at the first cap-proximal nucleotide (^{m7}Gppp^{m6}A_mN) is the second most frequently found modification in mRNA and occurs in conjunction with Cap 1 (and potentially Cap 2).

The immunogenic role of mRNA caps requires elucidation. Viral attenuation occurs after deleting methyltransferases that RNA viruses encode to convert Cap 0 to Cap 1. IFITs bind Cap 0 and activate antiviral translational repression. Thus, Cap 1 (and possibly Cap 2) marks endogenous mRNAs as "self" RNAs. The role of Cap 2 and Cap 1 (^{m6}A) is poorly understood because such capped mRNAs have not been produced synthetically at scale. In a recent study, Cap 1 (m6A) caps may increase stability and translation while decreasing de-capping of mRNA (Mauer et *a*l., Nature 2017, 541, 371-375).

Traditional co-transcriptional capping utilizes ARCA (Anti-Reverse Cap Analog) to produce immunogenic Cap 0 with poor capping (~70%) and low yield. Post-transcriptional enzymatic capping to produce Cap 0 or Cap 1 is hindered by highly structured 5' ends, requires further purification and is expensive. Methods to produce Cap 2 mRNAs have not been commercially available. We developed CleanCap[®], a novel cotranscriptional capping method to yield Cap 0, Cap 1, Cap 2, Cap 1 (^{m6}A) or unnatural caps (Figure). Capping with CleanCap is reproducibly efficient (90-99%), less expensive than enzymatic capping and is done in a "onepot" reaction without additional purification. In addition, CleanCap co-transcription method yields higher amount of capped mRNA than other methods including ARCA and enzymatic. Our studies in a THP-1 Dual monocyte cell line indicate that these various CleanCap mRNAs exhibit altered expression and immunogenicity. Further *in vivo* studies to characterize these mRNAs are ongoing.

Innate Immune Sensors Recognize mRNA

Transfection of cells with unmodified RNAs can lead to cell death due to activation of innate immune pathways Toll-like receptors 3, 7, & 8 recognize different RNA forms » Found in endosomes where some viruses enter cells

Cytosolic sensors

- » Protein Kinase R (PKR): dsRNA
- » MDA5: long dsRNA
- » IFITs: unmethylated cap structures » RIG-I: 5'-triphosphate

Background: Why mRNA Therapeutics?

- ▶ mRNA is a popular new tool for gene expression » Does not have a risk of insertional mutagenesis » Can transfect difficult cells such as non-dividing cells » Is transient
- Applications
- » Genome editing (Transposons, Cre, ZFNs, TALENs and CRISPR/Cas9) » Gene replacement » Vaccines
- Limitations
- » Innate immune response to unmodified mRNA
- Solutions
- » Proper capping
- » Chemical modification of mRNA can prevent innate immune stimulation » Removal of dsRNA

Figure 1: Cap 0, Cap 1 and Cap 2 Structures of 5'-Ends of mRNAs

Eukaryotic mRNAs have a Cap 1 or Cap 2 structure.

Sensing of proper cap structure is thought to be involved in self/non-self RNA recognition.

Cap 0: \mathbb{R}^1 and \mathbb{R}^2 =H Cap 1: $R^1 = CH_3$; $R^2 = H$ Cap 2: $R^1 = CH_3$; $R^2 = CH_3$ $B^{1-3} = A, C, G, U \text{ or } {}^{m6}A (B^1)$



Co-transcriptional capping with CleanCap (Cap 1) helps evade an immune response

» Splicing » Turnover » De-capping 2.7E+007 2.2E+007 1.6E+007 1.1E+007 5.4E+006 0.0E+000



¹TriLink BioTechnologies LLC, San Diego, CA 92121, USA ²Synthetic Genomics, Inc., La Jolla, CA 92037, USA







Capping Method	Capping Efficiency	Cost/mg Capped RNA	Cap 2 Achievat
Enzymatic	Variable	High	No
ARCA	~70 %	Moderate	No
CleanCap	~90-99%	Low	Yes

Figure 9: Protein Expression for Cap 0, Cap 1, Cap 2 or Cap 1 (^{m6}A) HPLC Purified Luciferase mRNAs

a) At 6 and 24 hours Luciferase Protein was Measured in

Cap 0 is inferior to Cap 1, Cap 2, and Cap 1 (^{m6}A) at 6 and 24 hours. Note: ARCA also yields ^{m6}A_m mRNAs are significantly superior to Cap 1 and Cap 2 mRNAs at 6 and 24 hours.

b) In Vivo Time-course of Luciferase Expression in Liver-

Not HPLC purified 5-MoU is showed similar activity to HPLC purified WT mRNA

c) In Vivo Time-Course of Luciferase Expression in Liver-

Prilink BIOTECHNOLOGIES part of Maravai LifeSciences

Figure 10: In Vitro De-Capping with Dcp2 is Decreased with ^{m6}A Capped RNAs. Identity of First and Second Cap-Proximal Influences De-capping

³²P labeled capped oligonucleotides with different cap forms and 5' sequences were decapped in vitro with purified capping enzymes

Conclusions for Cap 1^{m7}Gppp^{m6}A_mG is de-capped more slowly than Cap 1^{m7}GpppA_mG » The identity of the second cap proximal nucleotide ^{m7}Gppp^{m6}A_mN influences the rate of Dcp2 mediated de-capping

» ^{m6}AG displays the lowest de-capping rate for all capped forms

Figure 11: Nudt12 and DXO Selectively De-Cap Cap 0 but Not Cap 1 or Cap 2 RNAs

Data courtesy of Samie Jaffrey (Cornell) and Mike Kiledjian (Rutgers)

Conclusions

» CleanCap is a novel co-transcriptional capping method

- » Very high and consistent capping efficiencies obtained with CleanCap » CleanCap is an attractive, cost effective alternative to enzymatic or ARCA capping of mRNA
- » CleanCap allows novel cap forms that were not previously accessible such as Cap 2 and Cap 1 (m6 A)
- » Cap 1 and Cap 1 (^{m6}A) Cap RNAs are more active than Cap 0 RNAs *in vivo* » Cap 1 (^{m6}A) Cap alters activity *in vivo* and may extend persistence of ^{m6}A_m capped RNAs
- » ^{m6}A_mG RNAs are de-capped more slowly than A_mG RNAs
- » The identity of the second cap proximal nucleotide (${}^{m7}\text{Gppp}{}^{m6}\text{A}_{m}\text{N}$) influences the rate of Dcp2 mediated de-capping
- » ^{m6}A_mG displays the lowest de-capping rate for all capped forms
- » Nudt12 and DXO selectively de-cap Cap 0 but not Cap 1 or Cap 2 RNAs

Contact

Anton McCaffrey amccaffrey@trilinkbiotech.com

The Modified Nucleic Acid Experts® www.trilinkbiotech.com