

CleanAmp™ PCR 2X Master Mix

Catalog # L-5101

L-5101-100 (100 reactions)

L-5101-BK (Bulk amount)

CleanAmp™ PCR 2X Master Mix is an optimized, ready-to-use mix of CleanAmp™ dNTPs and *Taq* DNA Polymerase in reaction buffer for standard PCR. Simply add primers, template DNA and water.

QC Analysis

Functional Assay; Pass

Tested in standardized PCR assay for efficiency and specificity.

Handling & Use

Store at -20°C

Stable to 15 freeze-thaw cycles. Exposure to ambient temperatures during shipping does not adversely affect product performance.

Product released by Quality Assurance

CleanAmp™ Products: Patent Pending | RESEARCH LICENSE AGREEMENT

PURCHASE AND/OR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. CleanAmp™ products and/or technology incorporated therein (hereinafter "Products") and their use may be covered by one or more patents or pending Patent Applications. If you do not agree to use the Products you have purchased pursuant to the terms and conditions set out in this Research License Agreement, please contact TriLink BioTechnologies, Inc. within ten days of receipt to return the unused and unopened Products for a full refund; provided, however, that custom-made Products may not be returned for a refund.

1. **Research Use.** The purchase of Products conveys to the buyer the non-transferable right to use the purchased amount of the Products in internal research conducted by the buyer, whether the buyer is an academic, non-profit, or for-profit entity. Buyer agrees that it will not sell or otherwise transfer Products, or any components or derivatives thereof, to any third party. Notwithstanding the foregoing, materials made through use of the Products may be transferred by you to your legal affiliates or bona fide third party contractors performing paid work on your behalf, provided the use by such third party contractors is limited to performance of work for you and such work is performed subject to the terms of this Research License Agreement.

2. **Commercial Use.** Buyer also agrees that it will not sell, transfer, or otherwise use Products, or any components or derivatives thereof, for any commercial purposes. A license is required for any commercial use of Products, components, or derivatives thereof, regardless of the academic or non-profit status of the using entity. Information about commercial licenses for Products may be obtained by emailing info@trilinkbiotech.com or by writing to TriLink BioTechnologies, Inc., 9955 Mesa Rim Road San Diego, CA 92121. You may not use the Products to support the filing of a patent application in any country in the world that contains claims directed to the Products or uses thereof without the express approval of TriLink BioTechnologies, Inc.

3. **Attribution.** Users of the Products will expressly refer to the provision of the product in their published and unpublished works by explicitly identifying the Products purchased and stating that the Products were "purchased from TriLink BioTechnologies, Inc. (www.trilinkbiotech.com)."

4. **WARRANTY.** The Products are provided without warranty of merchantability or fitness for a particular purpose or any other warranty, express or implied, and without any representation or warranty that the use or supply of the Products will not infringe any patent, copyright, trademark or other right. TriLink BioTechnologies, Inc. does not recommend to its end users any particular application, methodology and/or protocol for the use of the Products. Depending on your particular use of the Products, it may be necessary to obtain a separate license or licenses from one or more third parties. TriLink BioTechnologies, Inc. and its employees and agents shall not be held liable for your use of the Products transferred to you. You agree to hold TriLink BioTechnologies, Inc. and its employees and agents harmless for any loss, claim, damage or liability, of whatsoever kind or nature, which may arise from acceptance, use, handling or storage of the Products by you.

5. **LIMITATION OF LIABILITY.** In no event shall buyer be entitled to recover from TriLink BioTechnologies, Inc. any special, indirect, incidental, consequential, or punitive damages in connection with this agreement, buyer's use of the Products, or the license granted hereunder.

6. Upon receipt of Products, buyer shall use its expertise and facilities in strict accordance with all applicable local, state and federal laws, regulations and guidelines. Buyer understands that the Products may have biological and/or chemical properties that are unpredictable and unknown at the time of transfer, that they are to be used with caution and prudence, and that they will not be used for testing in, or treatment of, humans.

7. Your right to have and use the Products will terminate immediately if you fail to comply with these terms and conditions of this agreement. You shall, upon such termination of your rights, destroy all Products, or any components or derivatives thereof, and notify TriLink BioTechnologies, Inc. of such in writing.

8. This Research License Agreement sets forth the complete and entire agreement of the parties with respect to the subject matter hereof and supersedes and terminates all prior agreements and understandings between the parties. No subsequent amendment or addition to this Research License Agreement shall be binding upon the parties unless reduced to writing and signed by the respective authorized officers of the parties. This Research License Agreement shall not be assigned or otherwise transferred by the buyer.

Protocols

Endpoint PCR (25 µL)

1. Thaw CleanAmp™ PCR 2X Master Mix, primers and DNA template and place on ice.
Note: Do not vortex CleanAmp™ PCR 2X Master Mix. Mix thoroughly by pipetting up and down and collect by pulse centrifugation.
2. Prepare a reaction mixture containing all components except for the DNA template. Add CleanAmp™ PCR 2X Master Mix, primers and sterile de-ionized water as shown in Table 1 into thin-walled PCR tubes. Keep on ice.
3. Mix the reaction mixture gently to protect the enzyme, by pipetting up and down. Do not vortex. Pulse spin if necessary.
4. Add the appropriate volume of template DNA to reach a reaction volume of 25 µL.
5. Pulse spin to remove bubbles and collect reaction solution at bottom of PCR tube.
6. Place the tubes into a thermal cycler with a heated lid and perform the appropriate cycling conditions for standard thermal cycling:
95°C for 0-5 min
[95°C for 10-40 sec; 48-60°C¹ for 1-30 sec; 72°C for 0.5-2 min²]
30-40 cycles, 72°C for 10 min
7. Analyze an aliquot of the completed reaction by agarose gel electrophoresis.

Table 1

Component	Final Concentration (25 µL reaction)	Volume per reaction
CleanAmp™ PCR 2X Master Mix	1X	12.5 µL
Forward/Reverse Primer	50-500 nM	Variable
DNA Template	Variable	Variable
Sterile De-ionized Water	Up to 25 µL	Up to 25 µL
Total Volume (µL)	25 µL	25 µL

¹ The annealing temperature should be chosen for optimal PCR performance. Most primer design software recommends an annealing temperature. The annealing temperature can also be optimized experimentally by using a thermal cycler with gradient functionality or by performing sequential experiments in which the annealing temperature is varied.

² The extension time at 72°C is recommended to be 30-60 seconds per kb of target.

Real-Time PCR

CleanAmp™ PCR 2X Master Mix has been successfully adapted for real-time detection using intercalating dye and probe-based detection. Please refer to the instrument manufacturer for specific protocols.

Real-time PCR may be proprietary. No license is conveyed expressly or by implication to the purchaser by purchase of any TriLink BioTechnologies products.

See Reverse for Troubleshooting

Troubleshooting

Probable Cause	Suggestion(s)
No amplification product or low amplicon yield	
Insufficient activation of CleanAmp™ dNTPs during thermal cycling	Optimize the duration of the initial denaturation time to up to 10 minutes.
Thermal cycling protocol is not optimized	Increase extension time. Generally extension times should be 30-60 seconds per kb of target.
	Increase the number of thermal cycles in 5 cycle increments.
	Optimize annealing temperature.
Problem with reagents or reaction conditions	Prepare fresh reagents, including CleanAmp™ PCR 2X Master Mix and primers.
	Verify that template is good in quality and of sufficient quantity.
	Verify primer design to ensure adequate complementarity to the DNA target.
Non-specific product formation	
Excessive off-target primer extension	Titrate the concentration of the primers or template DNA.
Primer dimer formation	Reduce initial denaturation and denaturation times: 95°C for 0-5 min, [95°C for 10-20 sec, 48-60°C for 1-15 sec, 72°C for 0.5-2 min] 30-40 cycles, 72°C for 10 min Note: A zero initial denaturation time in primer/template systems prone to primer dimer formation may cause a slight delay in threshold cycle (Ct) or quantification cycle (Cq).
Mis-priming	Omit initial denaturation time and shorten annealing time: [95°C for 30 sec, 48-60°C for 1 sec, 72°C for 0.5-2 min] 30-40 cycles 72°C for 10 min

Contact

TriLink BioTechnologies, Inc.
Info@trilinkbiotech.com
9955 Mesa Rim Road | San Diego, CA 92121
P(858) 546.0004 :: (800) 863.6801 | F(858) 546.0020