

CleanAmp™ dUTP Set

Catalog # N-9508

N-9508-2 (2 μmoles of each component)

N-9508-10 (10 μmoles of each component)

Sodium Salt form | Store at -20°C or below

Catalog Number	Component	Concentration
N-9520	CleanAmp™ dATP	50 mM
N-9521	CleanAmp™ dCTP	50 mM
N-9522	CleanAmp™ dGTP	50 mM
N-9524	CleanAmp™ dUTP	50 mM

QC Analysis

Functional Assay; Pass
 Tested in standardized PCR assay

AX-HPLC*; Pass RP-HPLC*; Pass

³¹P NMR*; Pass Mass Spec*; Pass

*Components analyzed separately

Product released by Quality Assurance

CleanAmp™ Products: Patent Pending | RESEARCH LICENSE AGREEMENT

PURCHASE AND/OR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS. CleanAmp™ products and/or technology incorporated therein (hereinafter "Products") and their use may be covered by one or more patents or pending Patent Applications. If you do not agree to use the Products you have purchased pursuant to the terms and conditions set out in this Research License Agreement, please contact TriLink BioTechnologies, Inc. within ten days of receipt to return the unused and unopened Products for a full refund; provided, however, that custom-made Products may not be returned for a refund.

1. **Research Use.** The purchase of Products conveys to the buyer the non-transferrable right to use the purchased amount of the Products in internal research conducted by the buyer, whether the buyer is an academic, non-profit, or for-profit entity. Buyer agrees that it will not sell or otherwise transfer Products, or any components or derivatives thereof, to any third party. Notwithstanding the foregoing, materials made through use of the Products may be transferred by you to your legal affiliates or bona fide third party contractors performing paid work on your behalf, provided the use by such third party contractors is limited to performance of work for you and such work is performed subject to the terms of this Research License Agreement.

2. **Commercial Use.** Buyer also agrees that it will not sell, transfer, or otherwise use Products, or any components or derivatives thereof, for any commercial purposes. A license is required for any commercial use of Products, components, or derivatives thereof, regardless of the academic or non-profit status of the using entity. Information about commercial licenses for Products may be obtained by emailing info@trilinkbiotech.com or by writing to TriLink BioTechnologies, Inc., 9955 Mesa Rim Road San Diego, CA 92121. You may not use the Products to support the filing of a patent application in any country in the world that contains claims directed to the Products or uses thereof without the express approval of TriLink BioTechnologies, Inc.

3. **Attribution.** Users of the Products will expressly refer to the provision of the product in their published and unpublished works by explicitly identifying the Products purchased and stating that the Products were "purchased from TriLink BioTechnologies, Inc. (www.trilinkbiotech.com)".

4. **WARRANTY.** The Products are provided without warranty of merchantability or fitness for a particular purpose or any other warranty, express or implied, and without any representation or warranty that the use or supply of the Products will not infringe any patent, copyright, trademark or other right. TriLink BioTechnologies, Inc. does not recommend to its end users any particular application, methodology and/or protocol for the use of the Products. Depending on your particular use of the Products, it may be necessary to obtain a separate license or licenses from one or more third parties. TriLink BioTechnologies, Inc. and its employees and agents shall not be held liable for your use of the Products transferred to you. You agree to hold TriLink BioTechnologies, Inc. and its employees and agents harmless for any loss, claim, damage or liability, of whatsoever kind or nature, which may arise from acceptance, use, handling or storage of the Products by you.

5. **LIMITATION OF LIABILITY.** In no event shall buyer be entitled to recover from TriLink BioTechnologies, Inc. any special, indirect, incidental, consequential, or punitive damages in connection with this agreement, buyer's use of the Products, or the license granted hereunder.

6. Upon receipt of Products, buyer shall use its expertise and facilities in strict accordance with all applicable local, state and federal laws, regulations and guidelines. Buyer understands that the Products may have biological and/or chemical properties that are unpredictable and unknown at the time of transfer, that they are to be used with caution and prudence, and that they will not be used for testing in, or treatment of, humans.

7. Your right to have and use the Products will terminate immediately if you fail to comply with these terms and conditions of this agreement. You shall, upon such termination of your rights, destroy all Products, or any components or derivatives thereof, and notify TriLink BioTechnologies, Inc. of such in writing.

8. This Research License Agreement sets forth the complete and entire agreement of the parties with respect to the subject matter hereof and supersedes and terminates all prior agreements and understandings between the parties. No subsequent amendment or addition to this Research License Agreement shall be binding upon the parties unless reduced to writing and signed by the respective authorized officers of the parties. This Research License Agreement shall not be assigned or otherwise transferred by the buyer.

Product Overview

PCR amplification products can cause false positives in downstream detection. Employment of Uracil-N-glycosylase (UNG), which recognizes and removes uracil residues, is the most common method for carryover decontamination in PCR. This method also substitutes dUTP for dTTP during PCR, resulting in amplicons that are distinct from the DNA template.

TriLink's CleanAmp™ dUTP is a hot start version of this analog, preventing nucleotide incorporation until the heat activation step. The patented CleanAmp™ technology employs thermolabile protecting groups blocking DNA polymerase function until activation by heat or pH adjustment.

Handling

CleanAmp™ dNTPs are characteristically very similar to natural dNTPs and are stable in aqueous buffer at pH 8-10.5 for at least one year at -20°C. When stored incorrectly, the major point for degradation of both natural and CleanAmp™ dNTPs is the triphosphate chain. Much like natural dNTPs, extra care should be taken not to expose CleanAmp™ dNTPs to prolonged temperatures above 20°C. Exposure to ambient temperatures during shipping does not adversely affect product performance. We do not recommend exposure of the CleanAmp™ dNTP stock solution to more than 24 TOTAL hours at room temperature or subjecting CleanAmp™ dNTPs to more than 20 freeze-thaw cycles. If greater than ten freeze-thaw cycles are anticipated, distribute your stock solution accordingly into smaller aliquots.

Use

Optimal PCR results are achieved when CleanAmp™ dUTP is in three-fold excess of dATP, dCTP and dGTP. Use CleanAmp™ dUTP with standard dNTPs for clean targets and with CleanAmp™ dNTPs for problematic targets. CleanAmp™ dUTP has been evaluated in routine PCR. See reverse for standard thermal cycling protocol and dNTP mixture preparation instructions. Should the CleanAmp™ dNTP concentration or units of DNA polymerase be increased, the MgCl₂ concentration should also be increased. Our data shows good PCR performance over a primer concentration range of 0.05 μM to 0.5 μM, finding 0.2 μM to work well in most cases.

See Reverse for Protocols and Troubleshooting

Protocols for *Taq* DNA Polymerase

Standard Thermal Cycling: 25 µL Endpoint PCR

1. For all components except CleanAmp™ dNTPs and DNA polymerase, thaw the reaction components, vortex to mix, centrifuge briefly and store on ice.
2. Prepare CleanAmp™ dNTPs:
 - a. Thaw at room temperature or on ice.
 - b. Vortex and pulse centrifuge to thoroughly mix.
 - c. If necessary, remove an aliquot of the stock solution and dilute with water or buffer (pH 8-10.5) to desired working concentration.
3. Prepare a mastermix containing all components except for the DNA template sample. Add each of the components as shown in Table 2 (multiply amounts by the number of reactions needed) in a centrifuge tube, on ice.
4. Mix the mastermix gently to protect the enzyme, by pipetting up and down. (Do not Vortex.) Pulse spin if necessary.
5. Aliquot 20 µL of mastermix into each thin-walled PCR tube.
6. To each 20 µL aliquot of mastermix, add 5 µL of the appropriate template DNA for a final reaction volume of 25 µL.
7. Pulse spin PCR tubes to remove bubbles and collect reaction solution at bottom of tube.
8. Place the tubes into a thermal cycler with a heated lid and perform the appropriate cycling conditions for target amplification:
95 °C for 10 min
[95 °C for 40 sec; X °C¹ for 1-30 sec; 72 °C for 0.5-2 min]
35-40 cycles 72 °C for 10 min
9. Analyze an aliquot of the completed reaction by agarose gel electrophoresis.

Table 2

Component	Final Concentration (in a 25 µL reaction)	Volume for 1 reaction	Volume for 10 reactions
Forward/Reverse Primer	50-500 nM	Variable	Variable
Sterile de-ionized water	Up to 25 µL	Up to 25 µL	Up to 250 µL
MgCl ₂ (50 mM) ²	4.5 mM	2.25 µL	22.5 µL
PCR buffer 10X ³	1X	2.5 µL	25 µL
CleanAmp™ dUTP Mix ⁴	0.2 mM	1 µL	10 µL
<i>Taq</i> 5 U/µL	0.05 U/µL	0.25 µL	2.5 µL
Template	Variable	5 µL	50 µL
Total Volume (µL)	25 µL	25 µL	250 µL

¹ The annealing temperature should be chosen for optimal PCR performance. Most primer design software recommends an annealing temperature. The annealing temperature can also be optimized experimentally by using a thermal cycler with gradient functionality or by performing sequential experiments in which the annealing temperature is varied.

² For protocols utilizing alternate DNA polymerases, the MgCl₂ concentration can be optimized between 3.0 and 5.0 mM.

³ 10X PCR buffer for Invitrogen's *Taq* DNA polymerase (200 mM Tris-HCl, pH 8.4 at 25°C; 500 mM KCl)
CleanAmp™ dUTP should be used in a 3-fold excess over dATP, dCTP and dGTP. See CleanAmp™ dUTP Mix Preparation below for the recommended nucleotide composition.

CleanAmp™ dUTP Mix Preparation

Table 1

Component	Stock Concentration ⁶	Desired Concentration	Volume ⁶
dATP ⁵	50 mM	10 mM	4 µL
dCTP ⁵	50 mM	10 mM	4 µL
dGTP ⁵	50 mM	10 mM	4 µL
CleanAmp™ dUTP	50 mM	30 mM	12 µL
Sterile Deionized Water	-	-	16 µL
			40 µL

⁵ Standard or CleanAmp™ dNTPs can be used.

⁶ Typical stock concentrations shown. Adjust volumes accordingly, as needed.

Troubleshooting

- Should excess banding be seen; try using less DNA polymerase.
- The annealing temperature should be chosen for optimal PCR performance.
- For improved performance, the CleanAmp™ dNTP concentration can be increased. For every additional 0.2 mM concentration of CleanAmp™ dNTPs, add at least an additional 1.0 mM of MgCl₂.
- Use CleanAmp™ dUTP (N-9524) with standard dNTPs (N-2505) for clean targets or with CleanAmp™ dNTPs (N-9507) for more problematic targets.
- Titrate template or primer concentrations for better specificity.

Contact

TriLink BioTechnologies, Inc.
info@trilinkbiotech.com
9955 Mesa Rim Road | San Diego, CA 92121
P 858-546-0004 | 800-863-6801