

## Description

Adenosine-5'-triphosphate (ATP), cytidine-5'-triphosphate (CTP), guanosine-5'-triphosphate (GTP), and uridine-5'-triphosphate (UTP) are unmodified ribonucleoside triphosphates, serving as essential building blocks of RNA molecules. Each of them consists of an unmodified base, a ribose, and a 5' triphosphate group. They can be used with modified nucleotides and TriLink's CleanCap® analogs in RNA synthesis by *in vitro* transcription.

TriLink offers both research-use-only (RUO) and good manufacturing practice (GMP) rNTPs. TriLink's GMP rNTPs are manufactured in highly controlled environments with documented procedures, traceability of materials and processes, and strict quality control measures, providing exceptional consistency, purity, and safety.

TriLink's RUO and GMP rNTPs may be ordered using the following catalog numbers:

	RUO				GMP		
	10 µmol	25 µmol	100 µmol	bulk	1 mmol	10 mmol	bulk
<b>ATP</b>	N-1501-10	N-1501-25	N-1501-100	N-1501-bk	FN-1501-1mmol	FN-1501-10mmol	FN-1501-bk
<b>CTP</b>	N-1502-10	N-1502-25	N-1502-100	N-1502-bk	FN-1502-1mmol	FN-1502-10mmol	FN-1502-bk
<b>GTP</b>	N-1503-10	N-1503-25	N-1503-100	N-1503-bk	FN-1503-1mmol	FN-1503-10mmol	FN-1503-bk
<b>UTP</b>	N-1504-10	N-1504-25	N-1504-100	N-1504-bk	FN-1504-1mmol	FN-1504-10mmol	FN-1504-bk
<b>rNTP set</b>	N-1505-10	N-1505-25	N-1505-100	N-1505-bk	—	—	—

They are provided 100 mM in H<sub>2</sub>O, pH 7.5.

## Use & handling

Store at or below -20°C. Upon first use, prepare aliquots to minimize freeze-thaw cycles. Use only certified RNase-free reagents and consumables with proper RNase-free technique.

## QC analysis

RUO products

- pH
- Concentration: UV/Vis
- Purity: <sup>31</sup>P NMR, AX-HPLC
- Identity: <sup>1</sup>H NMR, mass spectrometry

GMP products (in addition RUO analyses)

- Characterization: appearance, residual chemicals
- Safety: endotoxin, bioburden
- Nuclease contamination: DNase and RNase detection

Product released by Quality Assurance. TriLink is certified ISO 9001:2015.

## Troubleshooting

For any questions or technical support around this product, please reach out to [support@trilinkbiotech.com](mailto:support@trilinkbiotech.com)

## Other products in this guide

CleanCap® Reagent M6, also known as CleanCap® m6AG 3' OMe, is a TriLink's patented co-transcriptional cap analog to produce mRNAs containing a naturally occurring Cap-1 structure with >95% capping efficiency. Cap-1 mRNAs show superior *in vivo* activity compared to Cap-0 mRNAs produced by legacy capping methods such as mCap or anti-reverse cap analog (ARCA).

mRNA constructs created with CleanCap® Reagent M6 may further increase protein expression relative to previous generations of cap analogs, such as CleanCap® AG or CleanCap® AG (3' OMe), or mRNAs produced by enzymatic capping strategies<sup>1</sup>.

Modified uridines such as N1-methylpseudouridine, 5-methoxyuridine, and pseudouridine can reduce immunogenic response and enhance translational efficiency of mRNAs. These properties can result in safer mRNA and increased protein expression.

## Related TriLink products

---

Adenosine-5'-Triphosphate, ATP (cat no. N-1501)  
Cytidine-5'-Triphosphate, CTP (cat no. N-1502)  
Guanosine-5'-Triphosphate, GTP (cat no. N-1503)  
Uridine-5'-Triphosphate, UTP (cat no. N-1504)  
rNTP Set: ATP, GTP, CTP, UTP (cat no. N-1505)  
rNTP Set: ATP, CTP, GTP, N1MePsUTP (cat no. N-1506)\*  
rNTP Set: ATP, CTP, GTP, 5moUTP (cat no. N-1507)  
rNTP Set: ATP, CTP, GTP, PsUTP (cat no. N-1508)

CleanCap® Reagent M6 (cat no. N-7453)  
CleanCap® Reagent AG (cat no. N-7113)  
CleanCap® Reagent AG (3' OMe) (cat no. N-7413)  
CleanCap® Reagent AU (cat no. N-7114)

N1-Methylpseudouridine-5'-Triphosphate (cat no. N-1081)\*  
5-Methoxyuridine-5'-Triphosphate (cat no. N-1093)  
Pseudouridine-5'-Triphosphate (cat no. N-1019)  
T7 RNA Polymerase (Alphazyme cat no. E057)  
Inorganic Pyrophosphatase (E. coli) (Alphazyme cat no. E051)  
Engineered RNase Inhibitor (Alphazyme cat no. E07)

CleanCap® M6 EGFP mRNA (N1MePsU) (cat no. L-8101)<sup>‡</sup>  
CleanCap® M6 FLuc mRNA (N1MePsU) (cat no. L-8102)  
CleanCap® M6 mCherry mRNA (N1MePsU) (cat no. L-8103)  
CleanCap® M6 Cas9 mRNA (N1MePsU) (cat no. L-8106)<sup>§</sup>  
CleanCap® M6 EPO mRNA (N1MePsU) (cat no. L-8109)

## Related TriLink services

---

TriLink offers RUO and GMP custom mRNA production services in addition to our catalog offerings. Visit our website [trilinkbiotech.com/mrna-cdm-services](https://trilinkbiotech.com/mrna-cdm-services) or contact us at [mrna-services@trilinkbiotech.com](mailto:mrna-services@trilinkbiotech.com) for more information.

<sup>1</sup>Final capping is dependent upon the CleanCap® Reagent, DNA template, and final mRNA sequence. Secondary structure due to RNA length and base composition can affect final capping efficiency, mRNA yield, and translation efficiency.

Products containing CleanCap® technology are for internal research use only. A license is required for commercial use of CleanCap® Reagent M6 and other CleanCap® products. For license restrictions and patent(s) information, refer to the Research License Agreement below.

\*Nothing herein shall be deemed to constitute the grant of any license or other right in respect of any third-party rights in relation to the use of any materials or products controlled by a third party, and no license, sublicense or other rights shall be created hereunder by implication, estoppel or otherwise in or to any such third-party rights.

<sup>‡</sup>CleanCap® M6 Cas9 mRNA (N1MePsU) and/or other products or technologies relating to the Cas System (collectively, the "Cas Products") are provided under a Limited License granted by the Broad Institute, the Massachusetts Institute of Technology, President and Fellows of Harvard College, University of Iowa, University of Tokyo and Rockefeller University to the Buyer of the Cas Products, conveying to the Limited Licensee the non-transferrable right to use the purchased amount of the Cas Products solely for internal, non-clinical research to be conducted by the Limited Licensee found in [trilinkbiotech.com/legal-notices](https://trilinkbiotech.com/legal-notices)

<sup>§</sup>EGFP products are sold with a Limited Use Label License under sublicense with Life Technologies, Inc. Please review Limited License at [trilinkbiotech.com/gfp-label-license](https://trilinkbiotech.com/gfp-label-license).

## Quick start protocol for *in vitro* transcription (IVT) with CleanCap M6


[trilinkbiotech.com/cleancap-reagent-m6.html](http://trilinkbiotech.com/cleancap-reagent-m6.html)

### RNase-free techniques

It is essential that all reagents be RNase-free. Use disposable RNase-free tubes and bottles. Surfaces and pipettes can be wiped down with RNaseZap to destroy RNases. When possible, use dedicated RNase-free pipettes. Avoid using pipettes that have been used for plasmid preparation using RNase A.

### Template design

Template design is an integral part of any transcription. CleanCap M6 is to be used with the initiating sequence 5' AG 3'. The figure below shows the correct T7 promoter sequence (underlined) and initiator sequence (italics) for CleanCap M6.



### Customer supplied materials

NOTE: All reagents must be RNase-free. Use recommended source or equivalent grade.

#### Required reagents

- DNA template
- Nucleoside-5'-Triphosphate (NTP) Set (TriLink cat no. N-1505)  
Also available individually for use with modified NTPs. See related Products for commonly used modified NTPs.
- T7 RNA polymerase (New England BioLabs cat no. M0251S)
- Inorganic Pyrophosphatase (yeast) (New England BioLabs cat no. M2403S)
- Murine RNase Inhibitor (New England BioLabs cat no. M0314S)
- 1 M Tris-HCl (pH 7.5), RNase Free (Invitrogen cat no. 15567-027)
- Dithiothreitol (DTT) (EMD Millipore cat no. 3860-5GM)
- Spermidine (Sigma-Aldrich cat no. 85558-1G)\*
- 1N HCl (Fisher Scientific cat no. SA48)
- 1 M Magnesium Chloride (Sigma-Aldrich cat no. 63069)
- UltraPure™ DNase/RNase-Free Distilled Water (Thermo Fisher Scientific cat no. 10977015)

#### Optional reagents

- RNaseZap™ RNase Decontamination Solution (Thermo Fisher Scientific cat no. AM9780)
- DNase I (RNase-free) (New England BioLabs cat no. M0303S)
- CaCl<sub>2</sub> (Calcium chloride solution, BioUltra, 1M) (Sigma-Aldrich cat no. 21115)
- RNeasy Kits (QIAGEN cat no. 74104 or 75144)

\*NOTE: Preparation of 5 M spermidine stock solution

Undiluted spermidine stock (~6.37 M) is a viscous liquid that is difficult to pipette. To facilitate easier handling, we recommend making a 5 M dilution by adding 296 µL of nuclease-free water to a new 1 g bottle of 85558-1G and vortexing for a minimum of 60 seconds. This 5 M stock solution can then be aliquoted and stored at -20°C for at least 6 months. If using a larger starting mass than 1 g, scale the amount of water added accordingly.

## 10x M6 IVT transcription buffer preparation

### Reagent preparation

Add reagents in the prescribed order to ensure efficient transcription and capping. Ensure each component is well-mixed and homogenous before use, which is especially important for viscous reagents like spermidine.

**Table 1: 10x M6 IVT transcription buffer**

Component	10x concentration	Example 1000 $\mu$ L 10x M6 transcription buffer
DNase/RNase-free water	N/A	185.8 $\mu$ L
Tris pH 7.5 (1M)	400 mM	400 $\mu$ L
HCl (1M)	150 mM	150 $\mu$ L
MgCl <sub>2</sub> (1M)	160 mM	160 $\mu$ L
DTT (1M)	100 mM	100 $\mu$ L
Spermidine (5M)	21.2 mM	4.24 $\mu$ L

NOTE: To enable the highest capping efficiency, this 10x M6 transcription buffer contains extra HCl outside the buffering range of Tris. Do not titrate the 10x M6 transcription buffer to achieve a target pH. Preparing this 10x M6 transcription buffer with HCl results in a low pH (<2); however, once assembled into full 1x IVT reaction with all components the pH is close to neutral (~6.8).

## Standard-yield protocol

### Reagent preparation

Add reagents in the prescribed order to ensure efficient transcription and capping. Ensure each component is homogenous before use. Store thawed enzymes on ice. These reaction conditions have been tested with templates up to 6 kb in length.

1. Add RNase-free water and NTPs to the reaction tube.
2. Add CleanCap M6 to the tube and vortex to mix. Spin briefly to collect liquid.
3. Add 10x M6 transcription buffer as prepared according to Table 1. Vortex. Spin briefly to collect liquid.
4. Add DNA template.
5. Add Murine RNase Inhibitor, Inorganic Pyrophosphatase, and T7 RNA Polymerase.
6. Mix well by flicking or inverting the tube 10 times and spin briefly to collect liquid.
7. Incubate at 37°C for 3 hours.

**Table 2: Standard-yield reaction components**

Component	Final concentration	100 $\mu$ L rxn
DNase/RNase-free water	Up to 100 $\mu$ L	Up to 100 $\mu$ L
ATP (100 mM)	5 mM	5 $\mu$ L
CTP <sup>2</sup> (100 mM)	5 mM	5 $\mu$ L
GTP (100 mM)	5 mM	5 $\mu$ L
UTP <sup>2</sup> (100 mM)	5 mM	5 $\mu$ L
CleanCap M6 (100 mM)	10 mM	10 $\mu$ L
10x M6 transcription buffer	1X	10 $\mu$ L
DNA template	50 or 25 $\mu$ g/mL <sup>3</sup>	5 $\mu$ g or 2.5 $\mu$ g <sup>3</sup>
Murine RNase Inhibitor (40 units/ $\mu$ L)	1 unit/ $\mu$ L	2.5 $\mu$ L
Inorganic Pyrophosphatase (0.1 units/ $\mu$ L)	0.002 units/ $\mu$ L	2 $\mu$ L
T7 RNA Polymerase (50 units/ $\mu$ L)	15 units/ $\mu$ L	30 $\mu$ L
<b>Total volume</b>		<b>100 <math>\mu</math>L</b>

<sup>2</sup> Modified NTP can be used in place of wild-type. If using a modified NTP, use at the same concentration as the replaced wild-type NTP.

<sup>3</sup> Final concentration of DNA template should be 50  $\mu$ g/mL for a plasmid template or 25  $\mu$ g/mL for a PCR template.

## High-yield pulse feed transcription protocol

1. Add RNase-free water and NTPs to the reaction tube.
2. Add CleanCap M6 to the tube and vortex to mix. Spin briefly to collect liquid.
3. Add 10x M6 transcription buffer as prepared according to Table 1. Vortex. Spin briefly to collect liquid.
4. Add DNA template.
5. Add murine RNase inhibitor, inorganic pyrophosphatase, and T7 RNA polymerase.
6. Mix well by flicking or inverting the tube 10 times and spin briefly to collect liquid.
7. Incubate at 37°C for 2 hours. We recommend keeping the 10x M6 transcription buffer and NTPs at 4°C after the initial reaction set up to avoid an additional freeze thaw.
8. Approximately 5 minutes before the 2-hour incubation point, prepare the Spike-In mix in a separate tube, making sure to include 10% overage to account for pipetting error.
9. Vortex to mix. Spin briefly to collect liquid.
10. Add the appropriate volume of Spike-In mix (35 µL of Spike-In per 100 µL initial volume IVT) to the original reaction tube. Pipette to mix. Cap the tube and mix well by flicking or inverting 10 times. Spin briefly to collect liquid.
11. Incubate pulse-fed reaction at 37°C for an additional 2 hours (4 hours total).

**Table 3: High-yield pulse feed transcription protocol**

Initial IVT mixture	Component	Initial IVT mixture concentration	Volume	
	DNase/RNase-free water	Up to 100 µL	Up to 100 µL	
	ATP (100 mM)	5 mM	5 µL	
	CTP <sup>2</sup> (100 mM)	5 mM	5 µL	
	GTP (100 mM)	5 mM	5 µL	
	UTP <sup>2</sup> (100 mM)	5 mM	5 µL	
	CleanCap M6 (100 mM)	10 mM	10 µL	
	10x M6 transcription buffer	1X	10 µL	
	DNA template	50 or 25 µg/mL <sup>3</sup>	5 µg or 2.5 µg <sup>3</sup>	
	Murine RNase Inhibitor (40 units/µL)	1 unit/µL	2.5 µL	
	Inorganic Pyrophosphatase (0.1 units/µL)	0.002 units/µL	2 µL	
	T7 RNA Polymerase (50 units/µL)	15 units/µL	30 µL	
	<b>Total volume</b>			<b>100 µL</b>
<b>Incubate at 37 °C, 2 hours</b>				
Spike-In mixture	Component	Concentration change after pulse feed	Volume	
	ATP (100 mM)	Additional 4 mM	5.94 µL	
	CTP <sup>2</sup> (100 mM)	Additional 4 mM	5.94 µL	
	GTP (100 mM)	Additional 4 mM	5.94 µL	
	UTP <sup>2</sup> (100 mM)	Additional 4 mM	5.94 µL	
	10x M6 transcription buffer	Additional 0.74X <sup>4</sup>	14.85 µL	
	<b>Spike-In + 10% overage for 100 µL IVT</b>			<b>38.61</b>
	<b>Volume of Spike-In to add to 100 µL IVT</b>			<b>35 µL</b>
<b>Final high-yield transcription volume</b>			<b>135 µL</b>	
<b>Incubate at 37°C, 2 hours</b>				

<sup>2</sup> Modified NTP can be used in place of wild-type. If using a modified NTP, use at the same concentration as the replaced wild-type NTP.

<sup>3</sup> Final concentration of DNA template should be 50 µg/mL for a plasmid template or 25 µg/mL for a PCR template

<sup>4</sup>The final reaction buffer will be 1.74X after addition of the Spike-In mix

## Post-transcriptional options

### Purifications

IVT reactions may be purified by any traditional methods such as lithium chloride precipitation or spin columns (for example, QIAGEN RNeasy mini or midi kit) for higher purity at small scales. A fixed reaction will typically result in 4-5 mg of RNA per mL of reaction using wildtype NTPs or N1-methylpseudoUTP whereas a pulse feed reaction typically produces 8-10 mg of RNA per starting mL of reaction at approximately 7.5 mg/mL concentration (crude).

### DNase treatment

DNase treatment may be used per vendor-recommended protocol following first purification (above) and followed with a second small-scale RNA clean-up OR by recommended one-pot IVT/DNase reaction below.

For best one-pot DNase results formulate the final IVT reaction in 2 mM CaCl<sub>2</sub> with 20 U/μg of template DNA used in IVT at a final dilution of 4.8x the ending IVT volume by water. Mix by inverting or gently flicking the tube and pulse spin to collect liquid. Incubate for 20 minutes at 37°C.

**Table 4: Example DNase reaction setup**

Component	Final concentration	Fixed 100 μL IVT	Pulse-fed 135 μL IVT
IVT reaction	N/A	100 μL	135 μL
CaCl <sub>2</sub> (200 mM)	2 mM	4.8 μL	6.5 μL
DNase (2 U/μL)	20 U/μg of template DNA	25 μL <sup>6</sup>	25 μL <sup>6</sup>
Water	Variable	350.2 μL	481.5 μL
<b>Total volume</b>	<b>4.8x IVT volume</b>	<b>480 μL</b>	<b>648 μL</b>

<sup>6</sup> Volumes shown here based on IVT reaction using 25 μg/mL DNA template

## CleanCap® products | Research license agreement

[trilinkbiotech.com/cleancap-research-license](http://trilinkbiotech.com/cleancap-research-license)

PURCHASE AND/OR USE OF THIS PRODUCT SHALL CONSTITUTE ACKNOWLEDGMENT AND ACCEPTANCE OF THESE TERMS AND CONDITIONS.

Products containing the CleanCap technology (hereinafter “Products”) and their use may be covered by one or more patents or pending Patent Applications. If Buyer does not agree to use the Products purchased pursuant to the terms and conditions set out in this Research License Agreement (“Agreement”), the Buyer should contact TriLink BioTechnologies, LLC within ten days of receipt to return the unused and unopened Products for a full refund; provided, however, that custom-made Products may not be returned for a refund.

- 1. Research Use.** The purchase of Products containing CleanCap conveys to the buyer a non-exclusive, non-transferrable right to use the purchased amount of Products in internal research conducted by the buyer, whether the buyer is an academic, non-profit, or for-profit entity. Buyer agrees that it will not sell or otherwise transfer Products, or any components or derivatives thereof, to any third party. Notwithstanding the foregoing, materials made through use of the Products may be transferred by Buyer to Buyer’s legal affiliates or bona fide third party contractors performing paid work on Buyer’s behalf, provided the use by such third party contractors is limited to performance of work for Buyer and such work is performed subject to the terms of this Agreement.
- 2. Commercial Use.** Buyer also agrees that it will not sell, transfer, or otherwise use Products, or any components or derivatives thereof, for any commercial purposes, including (a) any human, clinical or clinical trial use, including, without limitation, any administration into humans or any diagnostic or prognostic use; (b) any human germline modification, including modifying the DNA of human embryos or human reproductive cells; (c) any in vivo veterinary or livestock use; (d) the development, manufacture, distribution, importation, exportation, transportation, sale, offer for sale, marketing, promotion or other exploitation or use of the Patent Rights or a Product for or as a testing service, therapeutic or diagnostic for humans or animals; (e) Products that provide nutritional benefits and are regulated by a regulatory authority as a drug or biologic pursuant to Section 505 of the Federal Food, Drug, and Cosmetic Act of 1938, as amended, Section 351 of the Public Health Service Act of 1944, as amended, or any successor laws, or equivalent laws or regulations in jurisdictions outside the United States; (f) any agricultural use, including but not limited to the use or application in the cultivation, growth, manufacture, exportation, or production of any tobacco product; and (g) any use or application relating to gene drive, unless and until a license is obtained for such commercial use of Products, components, or derivatives thereof, regardless of the academic or non-profit status of the using entity. Information about commercial licenses for Products may be obtained by contacting TriLink BioTechnologies, LLC. Buyer may not use the Products to support the filing of a patent application that contains claims directed to the Products or uses thereof in any country in the world without the express approval of TriLink BioTechnologies, LLC.
- 3. Attribution.** Buyers of the Products will expressly refer to the provision of the Products in their published and unpublished works by explicitly identifying the Products purchased and stating that the Products were “purchased from TriLink BioTechnologies, LLC ([www.trilinkbiotech.com](http://www.trilinkbiotech.com)).”
- 4. Warranty.** The Products are provided without warranty of merchantability or fitness for a particular purpose or any other warranty, express or implied, and without any representation or warranty that the use or supply of the Products will not infringe any patent, copyright, trademark or other right. TriLink BioTechnologies, LLC does not recommend to its end users any particular application, methodology and/or protocol for the use of the Products. Depending on Buyer’s particular use of the Products, it may be necessary to obtain a separate license or licenses from one or more third parties.
- 5. Limitation of Liability.** TriLink BioTechnologies, LLC and its employees and agents shall not be held liable for your use of the Products transferred to you. Buyer agrees to hold TriLink BioTechnologies, LLC and its employees and agents harmless for any loss, claim, damage or liability, of whatsoever kind or nature, which may arise from acceptance, use, handling or storage of the Products by Buyer. In no event shall buyer be entitled to recover from TriLink BioTechnologies, LLC any special, indirect, incidental, consequential, or punitive damages in connection with this agreement, buyer’s use of the Products, or the license granted hereunder.
- 6. Regulation Compliance.** Upon receipt of Products, buyer shall use its expertise and facilities in strict compliance with all applicable local, state and federal laws, regulations and guidelines. Buyer understands that the Products may have biological and/or chemical properties that are unpredictable and unknown at the time of transfer, that they are to be used with caution and prudence, and that they will not to be used for testing in, or treatment of, humans.
- 7. Termination.** Your right to have and use the Products will terminate immediately if Buyer fails to comply with the terms and conditions of this Agreement. Upon such termination of rights, Buyer shall destroy all Products, or any components or derivatives thereof, and notify TriLink BioTechnologies, LLC of such in writing.
- 8. Miscellaneous.** This Agreement sets forth the complete and entire agreement of the Parties with respect to the subject matter hereof and supersedes and terminates all prior agreements and understandings between the Parties. No subsequent amendment or addition to this Agreement shall be binding upon the Parties unless reduced to writing and signed by the respective authorized officers of the Parties. This Agreement shall not be assigned or otherwise transferred by the buyer.