TRILINK BIOTECHNOLOGIES PRESENTS

Next-generation mRNA design – increasing mRNA potency with a new cap analog

Kate Broderick, PhD, Chief Innovation Officer, Maravai LifeSciences

May 2023



part of Maravai LifeSciences

## Agenda

technology

**02** Introduction of our newest cap analog - CleanCap<sup>®</sup> M6

**03** Economics of CleanCap<sup>®</sup>



#### In Vitro transcription



#### IVT is the process of DNA template-directed RNA synthesis in a tube (or bioreactor)

- RNA polymerase
- Transcription template
- Standard and/or modified NTPs
- Cap analogs (optional)
- Inorganic pyrophosphatase
- RNase inhibitor
- Reaction buffer





#### Mature mRNA structure





4

## mRNA cap structure functions

mRNA cap structures are involved in modulating:

5



Splicing

#### > Translation initiation

**De-capping/turnover** 

Self/non-self recognition



#### Endogenous cap structures are methylated





## Endogenous cap structures are methylated

![](_page_6_Figure_1.jpeg)

![](_page_6_Picture_2.jpeg)

B₁

ÒМе

0=P-0

Ô

 $B_2$ 

OH

RNA

### Endogenous cap structures are methylated

![](_page_7_Figure_1.jpeg)

![](_page_7_Picture_2.jpeg)

B₁

OMe

O =

B<sub>2</sub>

ÔMe

RNA

## Capping strategies for mRNA manufacturing

![](_page_8_Figure_1.jpeg)

![](_page_8_Picture_2.jpeg)

## **Co-transcriptional capping**

#### ARCA vs CleanCap Reagent

![](_page_9_Figure_2.jpeg)

![](_page_9_Picture_3.jpeg)

## Simplified mRNA manufacturing with CleanCap<sup>®</sup> technology

![](_page_10_Figure_1.jpeg)

![](_page_10_Picture_2.jpeg)

## Expanding the CleanCap<sup>®</sup> family | Evolution of the technology

![](_page_11_Figure_1.jpeg)

![](_page_11_Figure_2.jpeg)

![](_page_11_Picture_3.jpeg)

#### Impact of cap structure on mRNA quality attributes

Sample	5'Cap	IVT Reaction Yield	Integrity	Capping Efficiency	Relative dsRNA	Residual DNA	Residual Protein
	ARCA⁺* (Cap 0)	1.3 mg/mL	78.6%	80%	+++	1.54 ng/mg	< 2%
N1-methylpseudo- UTP modified Firefly luciferase (FLuc)	CleanCapAG (Cap 1)	4.7mg/mL	93.3%	97%	++	1.27 ng/mg	< 2%
	CleanCapAG 3' OMe (Cap 1)	4.4 mg/mL	94%	96%	++	1.22 ng/mg	< 2%

mRNA integrity by IP-RP-HPLC; Capping efficiency by LC-MS; dsRNA by J2 blot; Residual DNA by qPCR; Residual Protein by NanoOrange

![](_page_12_Picture_3.jpeg)

#### Capping Strategy Data

#### Luciferase mRNA Expression in Mice

![](_page_13_Figure_2.jpeg)

![](_page_13_Picture_3.jpeg)

#### Improved protein expression with modification on m7G - CleanCap® AG 3'OMe

![](_page_14_Figure_1.jpeg)

![](_page_14_Figure_2.jpeg)

Integrated luminescent signal from N=5 animals across all time points displayed as mean with standard deviation. Two-tailed student's t-test was performed for all comparisons. ARCA signal is significantly lower than all other mRNAs (\*\*\*p<0.0001, \*p<0.01). CleanCap Reagent AG (3' OMe) is significantly higher than CleanCap Reagent AG (\*\*p<0.001).

![](_page_14_Figure_4.jpeg)

![](_page_14_Picture_5.jpeg)

## CleanCap<sup>®</sup> Analogs in COVID-19 Vaccines - CleanCap AG 3'OMe

#### A prefusion SARS-CoV-2 spike RNA vaccine is highly immunogenic and prevents lung infection in non-human primates

<sup>10</sup> Annette B. Vogel, Isis Kanevsky, Ye Che, Kena A. Swanson, Alexander Muik, Mathias Vormehr, Lena M. Kranz,

#### Received: 12 May 2020 Revised: 1 June 2020 Accepted: 1 June 2020 DOI: 10.1002/amp2.10060

COMMENTARY

ADVANCED MANUFACTURING MANUFACTURING AMPROCESSING

## Rapid development and deployment of high-volume vaccines for pandemic response

Zoltán Kis $^1 © \mid$ Cleo Kontoravdi $^1 \mid$ Antu K. Dey $^2 \mid$ Robin Shattock $^3 \mid$ Nilay Shah $^1$ 

#### nature

Explore our content V Journal information V

nature > articles > article

Article | Published: 30 September 2020

## COVID-19 vaccine BNT162b1 elicits human antibody and $T_H 1$ T cell responses

Ugur Sahin 🖾, Alexander Muik, [...] Özlem Türeci

Nature 586, 594–599(2020) Cite this article 95k Accesses 21 Citations 732 Altmetric Metrics

#### Manufacturing of RNA

BNT162b1 incorporates a Good Manufacturing Practice (GMP)-grade mRNA drug substance that encodes the trimerized SARS-CoV-2 spike glycoprotein RBD antigen. The RNA is generated from a DNA template by in vitro transcription in the presence of 1methylpseudouridine-5'-triphosphate (m1\PTP; Thermo Fisher Scientific) instead of uridine-5'-triphosphate (UTP). Capping is performed co-transcriptionally using a trinucleotide cap 1 analogue ((m<sub>2</sub>-7<sup>3'-O</sup>)Gppp(m<sup>2'-O</sup>)ApG; TriLink). The antigen-encoding RNA contains sequence elements that increase RNA stability and translation efficiency in human dendritic cells<sup>13,14</sup>. The mRNA is formulated with lipids to obtain the RNA-LNP drug product. The vaccine was transported and supplied as a buffered-liquid solution for intramuscular injection and was stored at -80 °C.

![](_page_15_Picture_17.jpeg)

![](_page_15_Picture_18.jpeg)

## Agenda

3Y

**02** Introduction of our newest cap analog - CleanCap<sup>®</sup>M6

**03** Economics of CleanCap<sup>®</sup> technology

![](_page_16_Picture_4.jpeg)

### Synthesis, screening and testing of novel cap structures

#### Screening process

- Cap analog synthesis Library of over 65 unique cap modifications
  - Modifications to m7G, modifications to 5'ppp5', modifications to +1A, etc.
- Screening those caps in IVT for capping efficiency with T7 RNA Pol.

![](_page_17_Figure_5.jpeg)

#### *In vivo* testing Firefly luciferase mRNA with over 20 new cap analogs

- LNP formulation
- Tail vein injection (1mg/kg)
- Body weights tracked out to 96 hours post mRNA delivery
- Imaging at 3, 6, 9, 12, 24, 36, and/or 48 hours post mRNA delivery by whole body bioluminescence

![](_page_17_Figure_11.jpeg)

![](_page_17_Picture_12.jpeg)

![](_page_18_Picture_0.jpeg)

# Our newest cap analog

![](_page_18_Figure_2.jpeg)

#### **CleanCap<sup>®</sup> Reagent M6**

![](_page_18_Picture_4.jpeg)

#### **m6Am modified caps occur naturally** – Applying rational design principals

doi:10.1038/nature21022

#### Nature Vol. 257 September 18 1975

<sup>13</sup> Sight, Land Stardt, R., Borow, W. Bight, N. Grandt, S. Mark, H. L. and M. K. Startman, Sequences of Held.a cell mRNAs revealed a network of the startman start of the startman startman startman startman startman startman startman startman startman, J. Startman, J. Startman, J. Startman, Startman, J. Startman, need nucleoside and found that is in common to mRNA of hum and monocel in well as a least new virtual and the digestion of HeLa cell mRNAs have the attractures m(C4) type(C) Humps and m(C2) regressive for the end of the second second second second second second end of the second second second second second second end of the second second second second second end of the second the second seco Marwin, A., Weininski, A., and Buchanan, J. M., Proc. neuro. 62, 440–453 (1987).
 RAVS, S., Castino, A., and Grifbuchek, E. P., J. molec. Biol., 54, 85–102 (1970).
 Travers, A. A., Narava, 223, 1107–1110 (1969).

#### N6, O2'-dimethyladenosine a novel methylated ribonucleoside next to the

5' terminal of animal cell and virus mRNAs MODIFIED 5'-terminal structures of the type m7G(5')ppp(5')Nm have been identified in various viral<sup>1+4</sup> and cellular<sup>1+7</sup> mRNAs. The 5-terminal modification of mRNA may be carried out by specific guanyly-1 and methyltransferases<sup>2</sup> and the structure seems to be required for efficient *in vitro* translation<sup>40</sup>. Studies<sup>6</sup> graphed between Am and N<sup>5</sup>-iopenteryladenousine (A). The munual nucleosities was designated AM and terratively considered to be a derivative of aderosities since is compared on the start of the star

![](_page_19_Figure_5.jpeg)

The 5 energial alignmethodium of HeLz cell mRNA were include and ingenetic with  $p_i$  nucleas cells energisms that hydro-layer RNA to 5 "nucleotidism and is not inhibited by 2-O entrybribunc-clocations" and sensitivity in the sensitivity of RNA to 5 "nucleotidisms". The sitter consignated with pA on paper electroperson at pH 33, threft-waggeding that the pAPm along with a small amount of contaministing pAm was there treated with 1 N (CL at 100 °C) of theorem for fits asso-them treated with 1 N (CL at 100 °C) of theorem for fits asso-them treated with 1 N (CL at 100 °C) of theorem for fits asso-them treated with 1 N (CL at 100 °C) of theorem for fits asso-them treated with 1 N (CL at 100 °C) of theorem for the pAPm along with a small amount of contaministing pAm was them treated with 1 N (CL at 100 °C) of theorem for the contamine the sensitivity of the sensitivity of the arrhydra was drived by treatmagnetic of the Ads. Some appeared has a fit the 2 O-methylabolism derivative resultive fits at the 2 -domethylabolism derivative results from the ecourties operiment in the hylabolism the and interact or the sensitivity of the arrangement of the solidated from the ecourties operiment and results and treated initiality. Fig. 1. Paper chromatography of methylated ribonucleosides from 5'sterminal oligonucleotides of: a, HeLa cell mRNA; the , mouse L. cell mRNA; and c. a denovirus mRNA. Cells were labelled for 4 h with methyl-H methionine, and polyadenylated mRNA real soluted from the cytoplasm of polyo(U)-Sophrace chromatography as previously described<sup>1</sup>. Polyadenylated mRNA from the cytoplasm of admovtine-infected HeLa cells was user tom the cytoplasm of adenovirus-infected Hela cells was also selected by hybriditation to adenovirus DNN (RM and F Koczot, ungublished). After hydrolysis of Hela, and L. cell mRNAs with RNas 12, TI and A and adenovirus mRNA with 0.4 N KOH the 5'-terninal oligonacleotides were induced DEAE-seluluos chrom-shared. mRNAs with RNase T2, T1 and A and adenovirus mRNA with 0.4 N KOH the 5-terminal oligonacheotids were isolated by DEAE-cellulose chromatography as the -5 to -6 charge peak'. The destaled 5-terminal oligonacleotids were then digested to their constituent ribonucleosides by a combination of venom phophotestersa and alkalian phophatas. The digest was analysed by accending paper chromatography with solvent containing p-bulanel-concentrated NH\_H\_d0.(86:5:14). chromatographed as m<sup>6</sup>Ade (Fig. 3b). Further evidence that the modified base of A<sup>4</sup>m is m<sup>6</sup>Ade was obtained by labelling HeLa cell mRNA with 2,8-<sup>3</sup>Hadenosine. Combined P1 nuclease and alkaline phosphatase digestion followed by adsorption of the enzyme-resistant

© 1975 Nature Publishing Group

Um Ö Gm

Fraction no. Fig. 2. Thin-layer electrophoresis of A\*m. Methyl-#H-labelled A\*m was isolated by paper chromatography as in Fig. 1a and analysed by electrophoresis on a thin layer of cellulose in 1 M formic acid at 800 V for 70 min.

The 5'-terminal oligonucleotides of HeLa cell mRNA were

Origin

#### ARTICLE

#### Reversible methylation of m<sup>6</sup>A<sub>m</sub> in the 5' cap controls mRNA stability

Jan Mauer<sup>1</sup>, Xiaobing Luo<sup>2</sup>, Alexandre Blanjole<sup>3</sup>, Xinfu Jiao<sup>2</sup>, Anya V. Grozhik<sup>1</sup>, Deepak P. Patil<sup>1</sup>, Bastlan Linder<sup>1</sup>, Brian F. Pickering<sup>1</sup>, Jean-Jacques Vasseur<sup>1</sup>, Qluying Chen<sup>1</sup>, Steven S. Gross<sup>3</sup>, Ollvier Elemento<sup>6,5</sup>, Françoise Debart<sup>2</sup>, Megeritich Kligdin<sup>4</sup> & Samie A. Jaffrey<sup>1</sup>

Internal bases in mRNA can be subjected to modifications that influence the fate of mRNA in cells. One of the most Internal bases in mRNA can be subjected to modifications that infinence the fate of mRNA in cells. One of the most prevalent modifications have shown that this matched the subject of th

As a mapping concept in gase sequence meghation is that a learner set of modified in section is final interval with mRNA, and bit materiy modification constitute on optimizer of the materiy modification constitute on optimizer of the materiy modified in the materiy on the materiy of the materiy modified in the materiy on the materiy modified in the materiy on the materiy of the materiy modified in the materia of the materiy modified in the materiy on the materiy modified in the materiy on the materiy of the materiy modified in the materia of the materiy of the materia materia of the materia materia of the materia concept in some expression regulation is that a diverse set for ETO. Using our transcriptome-wide man of m<sup>6</sup>A we find that

interaction in the control of the second se that one additional methyl modification can be detected in up to on their position within mRNAs. We measured the change in m<sup>6</sup>A 30% of mRNA caps<sup>9</sup>. If the first nucleotide following the m<sup>2</sup>G cap is stoichiometry for each m<sup>4</sup>A peak mapped in the Fto-knochout rela-2'-O-methyladenosine (A<sub>40</sub>). It can be further methylated at the N<sup>4</sup>-tive to the wild-type transcriptome. We used a previously described witching the multiple store of the store of m<sup>4</sup>A containing 2-O-methylation is essentially always detected at the first nucleotide. Notably, the m<sup>4</sup>A stoichiometry in the *Pio*-knockout at the first nucleotide.

A or  $m^6A_m$  as the first nucleotide, but not A or  $m^6A_m$  is unknown. (Extended Data Fig. 1a). As err M (ref. 10). The function of  $m^A_{n,k}$  is unknown. Enter we show that the standed mRNA cap carrier dynamical meters of the standed mRNA cap carrier dynamical we recorridy showed that the antibodies used in these early mFA map-ring static band both cannih the band. As a rest of the physiological contrast of m (200 c on the readil) con-tend to m (300 k and (300 k m)). The meters of the mathematical static mathematical static mathematical static mathematical static that m  $m^A_{n,k}$  and in the mathematical static mathematical static mathematical static that m  $m^A_{n,k}$  and the mathematical static mathematical static mathematical static that m  $m^A_{n,k}$  and the mathematical static that m  $m^A_{n,k}$  and m  $m^A_{n,k}$  is the static that m  $m^A_{n,k}$  and m  $m^A_{n,k}$  is the static that m  $m^A_{n,k}$  and  $m^A_{n,k}$  is the static that m  $m^A_{n,k}$  and  $m^A_{n,k}$  is the static that  $m^A_{n,k}$  and  $m^A$ 

<sup>1</sup> Department of Themacology, Bell Conter Muchan, Carel Liverany, New York, New York 2005, U.X. <sup>1</sup> Department of Coll Instigue and Neuroscience. Targets Unceredy, Thematol Neuroscience and Annual Neuroscience and Annual Neuroscience and Annual Neuroscience and Neuroscience an

00 MONTH 2016 | VOL 000 | NATURE | 1 © 2016 Macmilian Publishers Limited, part of Springer Nature. All rights reserved

#### Molecular Cell Article

#### Differential m<sup>6</sup>A, m<sup>6</sup>A<sub>m</sub>, and m<sup>1</sup>A Demethylation Mediated by FTO in the Cell Nucleus and Cytoplasm

Jampbo Wal,<sup>1,2,7</sup> Fange, Lu,<sup>1,2,4</sup> Zhile, Lu,<sup>1,2,4</sup> Qill Fe, Ji,<sup>2,</sup> Xiaol B, Fe, Ji - P. Cody He, <sup>1,3</sup> Halling Shi, <sup>1,2</sup> Xiaolong Cui, <sup>1,2</sup> Pul Su,<sup>4</sup> Ame Rhungland, <sup>1,4</sup> Galarey Jul, <sup>1,4</sup> James Tech, <sup>1,4</sup> and Charal He<sup>1,2,4</sup>. <sup>1,4</sup> Topartneet of Chemistry, Department of Biochemistry and Molecular Biology, and Institute for Biophysical Dynamics, <sup>1,4</sup> Lei Vinestry, of Charage, 2028 Sat 37 Sates, <sup>1,4</sup> Charage, Lison

Howard Hughes Medical Institute, The University of Chicago, 929 East 57 Street, Chicago, IL 60637, USA

Tentile of Material Sciences, Weatlake Instituté for Advanced Stady, Weitlake University, 18 Silongahan Read, Hangshou 31064, Olina "Deastment of Systems Bloogh, Beckmen Bleard Instituté of Cling Hilps, Monray Instituté of Basic Medical Sciences, University for Solis, Tentilate of Medical Microbiology, Osio University Hospital, Rishnapstalei, Norway Instituté of Basic Medical Sciences, University of Osio, To Bos 1116 Blence, 1315 Dea, Norway

\*Synthetic and Functional Biomolecules Center, Beijing National Laboratory for Molecular Sciences, Key Laboratory of Bioorganic Chemistry and Molecular Engineering of Ministry of Education, College of Chemistry and Molecular Engineering, Peking University

Beijing 100871, China 7These authors contributed equally Lead Contact \*Correspondence: chuanhe@uchi

https://doi.org/10.1016/j.molcel.2018.08.011

#### SUMMARY

have been devoted to identify the substrate(s) and to understand FTO, the first RNA demethylase discovered, mediates the demethylation of internal N<sup>6</sup>-methyladeno-the demethylation of 3-methylthymine (3meT) in single-stranded sine (m<sup>6</sup>A) and N<sup>6</sup>, 2-O-dimethyladenosine (m<sup>6</sup>A<sub>m</sub>) at DNA (Gerken et al., 2007) and 3-methyluracii (3meU) in RNA (Jia sine (m\*A) and M\*, 2-0-dimethyladenosine (m\*A) at UNK (seaso if al.,2001) and 3-methylarad (seaso if al.,2001) and 3-methy that FTO binds multiple RNA species, including fication in mammalian mRNAs (Fulet al., 2014). Adjacent to the mRNA, snRNA, and tRNA, and can demethylate inter-nal m<sup>6</sup>A and cap m<sup>6</sup>A<sub>m</sub> in mRNA, internal m<sup>6</sup>A in U6 (Adams and Cony, 1975; Wei et al., 1976), with a portion of these RNA, internal and cap m<sup>6</sup>A<sub>m</sub> in snRNAs, and N<sup>1</sup>-meth-yladenosine (m<sup>4</sup>A) in tRNA. FTO-mediated demethy-1975a), deposited by a yet to be identified methyltansferase lation has a greater effect on the transcript levels of This modification was confirmed by transcriptome-wide mRNAs possessing internal m<sup>6</sup>A than the ones with m<sup>6</sup>A-seq and exists in considerably lower overall abundance cap  $m^6A_m$  in the tested cells. We also show that FTO can  $m^6A_m$  is known to be an *in vitro* substrate To an directly repress translation by catalyzing of FTO (Fu, 2012), with a recent study showing that m<sup>4</sup>A<sub>m</sub> stability ated RNA demethylation occurs to m<sup>6</sup>A and m<sup>6</sup>A<sub>m</sub> in mRNA and snRNA as well as m<sup>1</sup>A in tRNA INTRODUCTION

lizes mRNA by preventing DCP2-mediated decapping and microRNA-mediated mRNA degradation (Mauer et al., 2017). However, the functional relevance of m<sup>6</sup>A<sub>m</sub> removal by FTO has yet to be fully explored. Surprisingly, this study also suggested that internal mRNA m<sup>6</sup>As may not be relevant substrates of FTO (Mauer et al

Because of these intriguing phenotypes, extensive efforts

2017), despite reports of a range of biological processes affected Fat mass and obesity-associated protein (FTO) has been sug-gested to correlate with human obesity by genome-wide associ-FTO-mediated m<sup>6</sup>A demethylation is critical in the DNA UV ation study (GWAS) (Fawcett and Barroso, 2010). A genetic damage response, with the m<sup>6</sup>A methyltransferase complex variant of FTD has been shown to be associated with increased METL3/14 exhibiting the opposite function (Kiang et al., food intake (Cecil et al., 2008), whereas loss-of-function muta-2017); FTD plays noticeable roles in *Flaviviridae* family virus tions in FTO cause severe growth retardation (Boissel et al., infection by demethylating viral FNA m<sup>6</sup>A in host cells (Gokhal

www.nature.com/cr www.cell-research.com

CelPress

Cell Research 

#### LETTER TO THE EDITOR OPEN

Cap-specific, terminal N<sup>6</sup>-methylation by a mammalian m<sup>6</sup>Am methyltransferase

Cell Research (2019) 29:80-82; https://doi.org/10.1038/s41422-018-

suggesting that PCF1 is a specific methyltiansferase for the terminal m<sup>2</sup>Am. Encouraged by the in vico readit, we then expressed and pointed recombanity PCIP protein, and stated and the protein strength and the protein strength and substrates under in vitro conditions Supplementary informa-tion, Fig32.1 — Hughest activity of PCI was obtained with RM (Probe-1: beginning with a complete coge structure with coppen-ray and based detective activity was detected and an antiperformation of RM Probe-2 with an internal Am RM (RM and RM and R Dear Editor, Dynamic and reversible M<sup>6</sup>-methyladenosine (m<sup>6</sup>A) RNA methy-lation has been found to greatly inpact gene expression. leading to the field of epitranscriptomics.<sup>1</sup> Unlike m<sup>6</sup>A that is an internal modification, a terminal modification at mRNA cap in higher to the field or paramotificanci. Usake mTA that is an intermal modification, is streamed in model and on mTAL top in Spite (Tig Li Li). The first and competitions the second nucleotide after the "Amethydamic constraints the second nucleotide after the "Amethydamic constraints" and constraints the second nucleotide after the "Amethydamic constraints the second nucleotide after the "Amethydamic constraints" and and and and and and probability of the second nucleotide after the "Amethydamic constraints" and and and and and and and probability of the second nucleotide after the "Amethydamic constraints" and and and and and and and probability of the second nucleotide and the probability of the second nucleotide after the second nucleotide after the second nucleotide and the probability of the second the nucleotide the antifyth and the the discription of this modified the antifyth and the second probability of the second the probability of the second the the second nucleotide the antifyth and the second the second and modul the second probability of the seco

 tori the tagins with m<sup>2</sup> (cppe)Ar. We modified the purplication of the tagins with m<sup>2</sup> (cppe)Ar. We modified the purplication is provided in the tagent of tagent o route of cell lysates, based on the procedure originally reported

SERVICER NATURE

Received: 15 September 2018 Revised: 9 November 2018 Accepted: 10 November 2018 Published online: 28 November 2018

© The Author(s) 2018

that PCIE1 is a povel mammalian m<sup>6</sup>Am writer, which is specifi

![](_page_19_Picture_40.jpeg)

part of Maravai LifeSciences

## CleanCap<sup>®</sup> M6 – the newest cap analog in the CleanCap<sup>®</sup> portfolio

![](_page_20_Figure_1.jpeg)

part of Maravai LifeSciences

<sup>+</sup> Mauer, J., Luo, X., Blanjoie, A. et al. Reversible methylation of m6Am in the 5' cap controls mRNA stability. Nature 541, 371–375 (2017). https://doi.org/10.1038/nature21022

![](_page_20_Picture_3.jpeg)

## Growing the CleanCap<sup>®</sup> family | CleanCap<sup>®</sup> AG technology portfolio

![](_page_21_Figure_1.jpeg)

![](_page_21_Picture_2.jpeg)

## Maintenance of capping efficiency across CleanCap<sup>®</sup> family members

![](_page_22_Figure_1.jpeg)

Capping efficiency by LC-MS – all cap analogs produce >95% capping efficiency

W TriLink BIOTECHNOLOGIES part of Maravai LifeSciences

#### CleanCap<sup>®</sup> M6 analog quality attributes are maintained using different sequences and lengths

5' Cap	Constructs	Length	IVT Reaction Yield	Capping Efficiency	Relative dsRNA
	eGFP	~ 1kb	4.4 ± 0.1 mg/mL	97.9 ± 0.4%	++
CleanCap M6 N1-methylpseudo-UTP modified bases	FLuc	~2 kb	5.15 ± 0.1 mg/mL	99.6 ± 0.1%	+
	Cas9	~4.5 kb	5.1 ± 0.1 mg/mL	97.2 ± 1.2%	++

Capping efficiency by LC-MS; dsRNA by J2 blot

![](_page_23_Picture_3.jpeg)

## CleanCap<sup>®</sup> M6 analog significantly improves protein expression levels

#### M6 modification promotes higher protein translation (FLuc) than previously observed

![](_page_24_Figure_2.jpeg)

Performance of FLuc mRNA in an LNP-formulated, tail vein delivered mouse model. 1 mg/kg dose per group. Luciferase activity, as photons per second, is measured after luciferin injection. The difference between groups is cap analog structure. All other variables are controlled.

![](_page_24_Figure_4.jpeg)

![](_page_24_Picture_5.jpeg)

## CleanCap<sup>®</sup> M6 analog enables potential lower dosing scenarios

#### M6 modifications allows for lower mRNA dosing regimes to get similar protein expression levels

![](_page_25_Figure_2.jpeg)

All groups of significantly different, \*\*\* p < 0.001, one-way ANOVA Error bars are standard error of mean. n = 7/group

![](_page_25_Picture_4.jpeg)

### Capping strategy comparison for potent mRNA manufacturing

![](_page_26_Figure_1.jpeg)

![](_page_26_Figure_3.jpeg)

## mRNA quality attributes comparing CleanCap® technology to enzymatic capping

Sample	5' Capping	IVT Reaction Yield	Integrity	Capping Efficiency	Relative dsRNA	Residual DNA	Residual Protein
	Enzymatic (Cap 1)	4.4 mg/mL	92.8%	99%	++	1.27 ng/mg	< 2%
N1-methylpseudo- UTP modified Firefly luciferase (FLuc)	CleanCap AG (Cap 1)	4.7 mg/mL	93.3%	97%	++	1.12 ng/mg	< 2%
	CleanCap AG 3'OMe (Cap 1)	4.4 mg/mL	94%	96%	++	1.22 ng/mg	< 2%
	CleanCap m6 (Cap 1)	4.3 mg/mL	97%	97%	+	0.44 ng/mg	< 2%

mRNA integrity by IP-RP-HPLC; Capping efficiency by LC-MS; dsRNA by J2 blot; Residual DNA by qPCR; Residual Protein by NanoOrange

![](_page_27_Picture_3.jpeg)

## **CleanCap® M6 makes mRNA more potent** – superior protein expression compared to enzymatically capped mRNA

![](_page_28_Figure_1.jpeg)

![](_page_28_Figure_2.jpeg)

\*\*\* p < 0.001, two-tailed T test. Error bars are standard error of mean. n = 9/group

![](_page_28_Figure_4.jpeg)

\*\*\*

1.8E+12

![](_page_28_Picture_5.jpeg)

## Simplified mRNA manufacturing with CleanCap® technology

![](_page_29_Figure_1.jpeg)

TriLink BIOTECHNOLOGIES part of Maravai LifeSciences

## Available now

Standalone cap analog	Size	SKU	mRNA service	Modification	SKU
CleanCap Reagent M6	1 µmol	N-7453-1	CleanCap M6 capped custom mRNA	Unmodified Bases	L-7507
	5 μmol	N-7453-5		Custom Base Modified	L-7508
	10 µmol	N-7453-10		PsU	L-7509
	100 µmol	N-7453-100		N1-Me-PsU	L-7510

![](_page_30_Picture_2.jpeg)

### Clean Cap<sup>®</sup> M6 analog | The Clear Choice

![](_page_31_Picture_1.jpeg)

## **Higher protein expression**

- New cap structure can product 30%+
- Potential to increase potency of mRNA drug substance
- Lower doses results in higher manufacturing yield

#### **Capping efficiency**

![](_page_31_Picture_7.jpeg)

#### Provides category leading capping efficiency of >95%

- Increased IVT efficacy resulting in high manufacturing yield
- Demonstrates reduced immunogenicity compared to other cap analogs

#### Manufacturing ease

![](_page_31_Picture_12.jpeg)

#### Maintains the one-pot workflow benefit of CleanCap technology

- Simplified manufacturing process, decreasing process risk
- Lowers time, labor, and cost to manufacture

![](_page_31_Picture_16.jpeg)

## Agenda

01	History of capping technology
02	Introduction of novel cap analog – CleanCap®

03 Economics of CleanCap<sup>®</sup> technology

![](_page_32_Picture_3.jpeg)

M6

## Economics considerations of mRNA capping strategies

- Third party study on economics of mRNA manufacturing and cost drivers
- 30 interviews of Biopharmaceutical experts in North America and Europe
- Quantified the raw material, labor and opportunity costs and benefits of each strategy

![](_page_33_Picture_4.jpeg)

Reagent costs

![](_page_33_Picture_6.jpeg)

![](_page_33_Picture_7.jpeg)

Challenges at large scale manufacturing

![](_page_33_Picture_9.jpeg)

#### Reagent cost comparison

Estimated both enzymatic capping and ARCA technologies lead to higher total reagent cost per 1 g of produced mRNA, GMP-grade:

- ~\$248k for enzymatic,
- ~\$221k for ARCA
- ~\$215k for CleanCap<sup>®</sup> technology

![](_page_34_Figure_5.jpeg)

![](_page_34_Picture_6.jpeg)

#### Process development time

Process development time utilizing CleanCap co-transcriptional technology is estimated to be about ½ **the time required** to optimize enzymatic or ARCA capping

- Simplified manufacturing process
- Fewer failure points/lower risk profile

![](_page_35_Figure_4.jpeg)

![](_page_35_Figure_5.jpeg)

![](_page_35_Picture_6.jpeg)

## Total manufacturing cost comparison

When comparing the overall manufacturing costs of the three capping strategies, CleanCap<sup>®</sup> technology is expected to be **30%** less than enzymatic capping and 20% less than ARCA.

For 1 gram of GMP grade mRNA batch, this CleanCap technology showed an estimated savings of:

- ~\$135,000 compared to enzymatic capping
- ~\$110,000 compared to ARCA capping

![](_page_36_Figure_5.jpeg)

![](_page_36_Picture_6.jpeg)

#### Manufacturing Costs at Large Scale

**Costs are non-linear** in the range of 1 to 40 gram

At the greatest amount of production as predicted for Phase III studies,

ARCA and enzymatic capping are estimated to be 30% and 40% more expensive than CleanCap

![](_page_37_Figure_4.jpeg)

![](_page_37_Picture_5.jpeg)

## Conclusion slide – Benefits of using CleanCap® technology

![](_page_38_Picture_1.jpeg)

#### Reduce both reagent and labor costs

- Simplified and streamlined manufacturing process limited excess raw material consumption and labor costs
- Cost savings increase significantly as manufacturing processes scale into late clinical and commercial development

#### Cap1 structure – no extra steps

- CleanCap technology achieves Cap1 structure without the need for multiple enzymatic and purification processes
- Extra bioreactor and purification steps are a fixed cost per batch, even with column reuse

![](_page_38_Picture_8.jpeg)

#### Process development time is cut in half

- Simple processes decrease process risk and accelerate development timelines
- Preclinical velocity is key to bringing new therapies to market

![](_page_38_Picture_12.jpeg)

## Who we are | TriLink BioTechnologies

CleanCap® NTPs Determined to deliver innovative nucleic acid tools and **Products &** services to help you bring transformative nucleic acid Technology therapies from research to patients. **Technical expertise to support your programs:** Services Capping analogs, NTPs, and modified UTPs **Custom oligonucleotides** (~) **IVT** reaction Analytical pDNA mRNA **Custom chemistry** (~) Clinical GMP Plasmid manufacturing for mRNA (~) Discovery Expert manufacturing Support Discovery mRNA Services (>25 years, >975 customers served) (~) **Process Development, Analytics, Quality Systems** GMP mRNA Services (>100 batches to date) (~) Process Development, Quality Systems, and Analytical Services Technology Rational Innovation advances R&D

40 ©2023 Maravai LifeSciences. This document and its contents are confidential and proprietary and are not to be shared or redistributed without express consent.

part of Maravai LifeSciences

## **Our history** | Over 25 years of innovation in mRNA

![](_page_40_Figure_1.jpeg)

![](_page_40_Picture_2.jpeg)

## Facilities to meet the demand in mRNA products and services

![](_page_41_Picture_1.jpeg)

#### mRNA & related raw materials

- 118,000 ft<sup>2</sup> (10,963 m<sup>2</sup>)
- mRNA services for development and early clinical
- CleanCap<sup>®</sup> reagents and NTP innovation and scale up

![](_page_41_Picture_6.jpeg)

#### Nucleic acid production

- 32,000 ft<sup>2</sup> (2,973 m<sup>2</sup>)
- CleanCap Reagents and NTP manufacturing
- GMP raw materials for clinical and commercial use
- BARDA Award

![](_page_41_Picture_12.jpeg)

#### Late phase clinical mRNA manufacturing

- 32,000 ft<sup>2</sup> (2,973 m<sup>2</sup>)
- mRNA development and manufacturing space
- Phase 1 clinical and beyond

![](_page_41_Picture_17.jpeg)

## Acknowledgments

#### M6 Discovery and Synthesis

Chunping Xu Chanfeng Zhao Ilya Ilichev **Russell Cousins** Alexandre Lebedev Mike Houston

#### M6 Process Development

Evan Myers **Tran Trinh Rachel Schubert** Maya Onuma **Mitchell Sugii** 

#### M6 Product and Project Management

Matthew Poling Nina Schulte Jason Downing Cory Smith

Natalie Hauser

#### M6 Screening and IVT Optimization

Jordana Henderson Andrew Ujita Anthony Truong Coleen Vo Deidre Broxterman Farinaz Rezvani Nona Abolhassani **Stephanie Ramos** 

#### M6 Manufacturing and Process Engineering

Sarah Walker Holly Howarth **Kevin Huang** Jordan Eaton **Daniel Reitzler** Khatrin Rico **Ruben Gonzalez**  Paul Cativiela Amanda Maiello **Rick Tham** Vanessa Hernadez M6 Quality Control

Inna Koukhareva **Amy Kuzminov Brett Turney Michael Sierra** Mark Gianino Ideen Eslamian **Pascal Sodemann David Martinez** 

M6 Marketing

**Doreen Pippen** Jennifer Gelman **Joey Tarantino** 

![](_page_42_Picture_17.jpeg)

#### Contact us

To order CleanCap<sup>®</sup> M6 analog

sales@trilinkbiotech.com

CleanCap<sup>®</sup> Reagent M6 | TriLink BioTechnologies

## Thank you

![](_page_44_Picture_1.jpeg)

part of Maravai LifeSciences