

Key Critical Quality Attributes for the cGMP Production of Therapeutic Messenger RNA

Jessica Madigan, Craig Dobbs, Krist T. Azizian, Dongwon Shin, Sabrina Shore, Jordana M. Henderson, Alexandre Lebedev, Anton P. McCaffrey and Richard I. Hogrefe

TriLink BioTechnologies
San Diego, CA 92121



Abstract

Recently, there has been significant interest in the use of messenger RNA (mRNA) as an *ex vivo* or *in vivo* therapeutic. Since mRNA is expressed in the cytoplasm it may be particularly useful for improving gene expression in difficult to transfect non-dividing cells. In contrast to plasmid or viral vectors, there is no risk of insertional mutagenesis or subsequent oncogenesis upon mRNA transfection and the transient nature of mRNA expression is desirable for genome editing (CRISPR/Cas9, ZFNs and TALENs) and vaccines. In each case, the goal is to produce a synthetic RNA that mimics a natural mRNA.

In anticipation of clinical testing, it is essential to use sequence designs, raw materials and manufacturing processes that are scalable. Early on, critical decisions must be made about 1) transcription template designs, 2) capping methodologies, 3) sequence optimization, 4) chemical modifications to evade innate immune responses, 5) polyadenylation methods and 6) purification methods.

For optimal expression in cells or target organs, transfected mRNAs must avoid detection by pattern recognition receptors (PRRs) that evolved to sense improperly capped RNAs and double stranded RNA. PRR activation leads to cytokine production, translational arrest and cell toxicity or death. Mammalian mRNAs are modified post-translationally to contain modified nucleotides (e.g. pseudouridine and 5-methylcytidine). These modifications can reduce activation of PRRs and allow maximal translation.

During RNA capping, Cap 0 (m⁷GpppN) is formed as an intermediate. Methylation of the 2' position of the first nucleotide forms Cap 1 (found in 100% of transcripts) and Cap 2 (found in 50% of transcripts) of endogenous mRNAs. mRNAs generated with commercially available cap analogs (ARCA) contain Cap 0 structures that can be immunogenic. Recombinant enzymes used to generate Cap 1 mRNA are expensive, do not always go to completion and the RNA must be purified prior to capping. A novel co-transcriptional capping method that yields Cap 1 or Cap 2 with high efficiency and lower costs in a 'one pot' reaction will be discussed.

First generation mRNAs were modified with pseudouridine or 5-methylcytidine/pseudouridine and had Cap 0 structures. Data from improved second generation Cap 1 mRNAs containing a combination of sequence engineering and chemical modifications will be presented. We will provide a broad roadmap for the application of these principles to the design and manufacturing of novel mRNA therapeutics.

Function of mRNA Cap Structures


Eukaryotic mRNAs have a Cap 1 or Cap 2 structure.

Sensing of proper cap structure is thought to be involved in self/non-self RNA recognition.

Cap structure influences activation of PRRs

» RIG-I is activated by Cap 0 RNAs but not Cap 1 mRNAs (PMID: 18426922 and 20457754)

» IFT1 binds Cap 0 RNAs more tightly than Cap 1 mRNAs (PMID: 24371270)



Co-transcriptional capping with CleanCap™ (Cap 1) helps evade an immune response

Capping Methods

- » Enzymatic capping characteristics
 - Utilizes Vaccinia virus capping enzyme and 2' O-methyltransferase enzymes
 - Can potentially achieve 100 % capping with a natural unmodified cap structure
 - Very costly at scale with significant batch- to-batch enzymatic
 - Accessibility of 5' end important for efficient capping

Co-transcriptional capping

- » ARCA capping characteristics
 - Capping fairly reproducible (60-80 % capping, 20-40% uncapped)
 - There is a significant amount of transcriptional "stuttering"
 - Cap contains an unnatural 3' O-methyl group
 - Transcript must start with G

- » CleanCap™ 3rd generation cap analog
 - Capping efficiencies typically 90-99%
 - Can produce Cap 0, Cap 1 or Cap 2
 - Yields a natural unmodified cap structure
 - Allows any base (A, C, G, or U) at 5' end
 - Allows ^mA or ^mA_n at 5' end
 - More affordable than enzymatic capping




Figure 2: CleanCap™, Enzymatic and ARCA, Capping Comparison





Figure 3: Template Design Cas9

- » Endotoxin free
- » Kanamycin resistance
- » Epitope tags may be undesirable
- » Is cGMP plasmid required?
- » Cell banking
- » Think about plasmid sourcing early



mRNA Purification

Purification is required to remove: salts, NTPs, cap analogs, proteins, truncated mRNA products, residual DNA and sometimes double stranded RNA

TriLink has developed two complementary chromatographic mRNA purification methods:


LC-Isolation

- » Removes residual proteins, DNA, and NTPs
- » Highly scalable


Reverse phase HPLC (RP-HPLC)

- » Decreases double-stranded RNA
- » Enriches for full length mRNA
- » Longer sequences more challenging to purify

Figure 4: In Vitro Translation and Cell Activity of Modified Luciferase mRNAs



dsRNA slot blot HPLC vs. non-HPLC




Slot blot analysis of dsRNA present in mRNA samples. Protocol adapted from Pub Med ID 23296926.

Purification Conclusions

- » For some applications, LC-isolation may be of sufficient quality especially if using uridine-depleted 5-methoxyuridine-substituted mRNAs. May not require RP-HPLC
- » For other applications, RP-HPLC purification may be required to remove double-stranded RNA
- » RP-HPLC purification reduces interferon induction in THP-1 cells in many cases
- » When selecting synthesis scale, account for purification yield

Figure 5: In Vitro Translation and Cell Activity of Modified Luciferase mRNAs



Contact

Jessica Madigan
jmadigan@trilinkbiotech.com

The Modified Nucleic Acid Experts®
www.trilinkbiotech.com

TIDES Oligonucleotide and Peptide Therapeutics
Boston, MA
05/18

Background: Why mRNA Therapeutics?

- » mRNA is a popular new tool for gene expression
 - » Does not have a risk of insertional mutagenesis
 - » Can transfect difficult cells such as non-dividing cells
 - » Is transient

Applications

- » Genome editing (Transposons, Cre, ZFNs, TALENs and CRISPR/Cas9)
- » Gene replacement
- » Vaccines

Limitations

- » Innate immune response to unmodified mRNA

Solutions

- » Proper capping
- » Chemical modification of mRNA can prevent innate immune stimulation
- » Removal of dsRNA

Innate immune sensors recognize mRNA

Transfection of cells with unmodified RNAs can lead to cell death due to activation of innate immune pathways


Toll-like receptors 3, 7 & 8 recognize different RNA forms

- » Found in endosomes where some viruses enter cells

Cytosolic sensors

- » Protein Kinase R (PKR): dsRNA
- » MDAS: long dsRNA
- » IFITs: unmethylated cap structures
- » RIG-I: 5'-triphosphate

Figure 1: mRNA Synthesis by In Vitro Transcription



Schematic representation of transcription and translation with N1-modified pseudouridine 5'-triphosphate derivatives in the presence of CleanCap™. *N1-modified-Ψ