

CleanCap™ Cas9 mRNA

L-7606-1000

Catalog # (Pack Size): L-7606-20 (20 ug)

L-7606-100 (100 ug)

L-7606-1000 (1 mg)

L-7606-BK (bulk amount)

Lot: WOTL35238

RNA Length: 4521 NT

Modifications: Capped (Cap 1) using CleanCap™ Technology.

Packaged As: 1.00mg/mL in 1 mM sodium citrate pH 6.4

Quality Control Data

QC Release Date: 09/24/2019

UV-Vis (A260/A280 Ratio) = 2.21

Gel Analysis of Final Product - Pass

Bioanalyzer - Pass

Recommended Storage: -40°C to -80°C

Avoid repeated freeze/thaw cycles.

Upon first use, it is recommended to aliquot sample into single use portions.

Approved By: _____

Quality Assurance

Date

al Aquino

24 Sept 2019
