



Product Data Sheet

MagnaLink Amino Magnetic Beads

Storage: Store at 4°C



Catalog Number: M-1005	Diameter: 2.8 µm
Lot Number:	Expiration:
Primary amine content: XX.XX nmol/mg beads quantified by reaction with Sulfo-S-4FB	Density: 1.9 mg/mL
Solids content:	Beads per ml: 8×10^8
Components: Uniform polymer encapsulated magnetic bead containing primary amino functional groups (-NH ₂).	Size: 1.0 mL Units: 10 mg

Product Description

A 2.8 micron uniform polymer encapsulated magnetic bead containing primary amine functional groups

Application

Amino functional groups can be activated with any water-soluble NHS-activated linker.

Comments

Not for internal or external use in humans.

