


Ara-GTP

N-1100-1



Aruanosine-5'-Triphosphate

Product Information

Catalog # (Pack Size): N-1100-1 (1 μ mole)

N-1100-5 (5 μ mole)

N-1100-10 (10 μ mole)

N-1100-BK (Bulk Amount)

Lot #: WOTL46421

Extinction Coefficient: 13,600 $\text{Lmol}^{-1}\text{cm}^{-1}$ @ 252 nm

Molecular Weight: 523.1 g/mole (free acid)

Molecular Formula: $\text{C}_{10}\text{H}_{16}\text{N}_5\text{O}_{14}\text{P}_3$ (free acid)

Final Salt Form: Li^+

Handling: Avoid repeated freeze/thaw cycles.
Upon first use, it is recommended to aliquot sample into single use portio

Packaged As: 100 mM in H_2O

Recommended Storage: -20°C or below

QC Analysis

Method of Analysis:

AX-HPLC

UV-Scan

Mass Spec

^{31}P NMR

^1H NMR

Conductivity

Results:

Passes

Passes

Passes

Passes

Passes

Passes

QA Approval

Approved By:

QA

06/10/2020

Date